



TARIQ MAHMOOD
● ● ● ● ● enterprises



DENTAL INSTRUMENTS CATALOGUE

www.tmentp.com



TABLE OF CONTENTS

Page 01-05

Scissors

01

Page 06-14

Tweezers, Forceps & Needle Holder

02

Page 15-33

Extracting Forceps

03

Page 34-73

Mirror Handles, Scalars & Knives

04

Page 74-83

Root Elevators

05

Page 84-110

Implantology Instruments & Mouth Gags

06



07

Dental Syringes

Page 111 - 115

08

Crown, Ligature & Laboratory Instruments

Page 116 - 129

09

Articulators

Page 130 - 131

10

Orthodontic Instruments

Page 132 - 149

11

Impression Trays

Page 150 - 154

12

Hollow Ware

Page 155 - 164


13

Surgical Burs

Page 165 - 170

Page 01-05

- SCISSORS
- GUM & OPERATING SCISSORS
- SCISSORS WITH TUNGSTEN CARBIDE INSERTS
- MICRO SCISSORS



TM-01-100
11.5 cm
TM-01-102
14.5 cm



TM-01-101
11.5 cm
TM-01-103
14.5 cm



TM-01-104
11.5 cm
TM-01-106
14.5 cm



TM-01-105
11.5 cm
TM-01-107
14.5 cm



TM-01-108
11.5 cm
TM-01-110
14.5 cm



TM-01-109
11.5 cm
TM-01-111
14.5 cm



TM-01-112



TM-01-113



Standard
Sharp / Sharp



Standard
Blunt / Sharp



Standard
Blunt / Blunt



Mayo
14.5 cm



TM-01-114



TM-01-115



TM-01-116



TM-01-117




Metzenbaum-Delicate
14.5 cm



Metzenbaum
14.5 cm



TM-01-118
Littauer
14 cm



Spencer
TM-01-119 9 cm
TM-01-120 13 cm

Suture Scissors




Bandage Scissors Lister
TM-01-121 9 cm
TM-01-122 11 cm
TM-01-123 14 cm
TM-01-124 18 cm




Gum & Operating Scissors



Iris
11.5 cm



Noyes / Castroviejo
12 cm



Goldman-Fox
13 cm



TM-02-135
Universal
18 cm



TM-02-136
Kelly
16 cm

Kelly
16 cm



TM-02-139
Goldman-Fox
13 cm

TM-02-140
La Grange
11.5 cm



Scissors, One Blade Serrated

- One blade is with fine serrations to avoid tissue or ligature slippage.



Iris
11.5 cm



Goldman-Fox
13 cm




TM-03-147
Goldman-Fox
13 cm



TM-03-149



TM-03-148



Kelly
16 cm




TM-03-150
Locklin
16 cm



TM-03-151
Locklin
16 cm

Super-Cut

- Cutting Blade horned to "razor" sharpness with opposing blade finished with fine serrations to avoid tissue slippage.



TM-03-152


TM-03-153



Iris
11.5 cm



TM-03-156
La Grange
11.5 cm



TM-03-154


TM-03-155



Goldman-Fox
13 cm



Scissors With Tungsten Carbide Inserts



TM-04-160

TM-04-161



Iris
11.5 cm




TM-04-162

TM-04-163



Metzenbaum
14.5 cm




TM-04-164


TM-04-165



Kelly
16 cm



TM-04-159
La Grange
11.5 cm




TM-04-166

TM-04-167



Goldman-Fox
13 cm



TM-04-168

TM-04-169




Strabismus
11.5 cm



TM-05-170
15 cm



TM-05-171
15 cm



TM-05-172
18 cm



TM-05-173
18 cm



TM-05-174
18 cm



Micro Scissor
Stainless Steel



Micro Scissor
Stainless Steel




Micro Scissor
Stainless Steel
for ligature, one blade serrated

Micro Scissors with Titanium Stainless Steel Tips




TM-05-180
18 cm




TM-05-176
18 cm



TM-05-177
18 cm



TM-05-178
15.5 cm



TM-05-179
15.5 cm



Page 06–14

- TWEEZERS
- DRESSING & DISSECTING FORCEPS
- HAEMOSTATIC FORCEPS
- NEEDLE HOLDERS
- NEEDLE HOLDERS WITH TUNGSTEN CARBIDE INSERTS
- MICRO NEEDLE HOLDERS & FORCEPS



TM-06-181
College
15 cm

TM-06-182
College
15 cm

TM-06-183
College-Locking
15 cm

TM-06-184
Flagg
16 cm

TM-06-185
Meriam
16 cm

TM-06-186
College
15 cm

TM-06-187
Flagg
16 cm

TM-06-188
Meriam
16 cm



TM-06-189
Perry
13 cm

TM-06-190
Perry
13 cm



TM-06-191



TM-06-192



TM-06-193



TM-06-194




Gerald
18 cm



Gerald
18 cm




Dissecting & Tissue Forceps






Standard
Dressing Forceps



Standard
Dissecting Forceps



Tissue Forceps

Splinter Forceps


Towel Forceps



Adson-Brown
TM-08-237
12 cm



Stieglitz



Backhaus
TM-08-241 9 cm
TM-08-242 11 cm
TM-08-243 13 cm




Jones
TM-08-244
9 cm






Schaedel
TM-08-245
9 cm



Haemostatic Forceps







Needle Holders




TM-11-276
Derf
12.5 cm

TM-11-277
Olsen-Hegar
14 cm

TM-11-278
Baumgartner
14.5 cm

TM-11-279
Crilewood
15 cm

Mayo-Hegar
TM-11-280 14 cm
TM-11-281 15 cm
TM-11-282 16 cm
TM-11-283 18 cm
TM-11-284 20 cm



Mathieu
TM-11-285 14 cm
TM-11-286 17 cm
TM-11-287 20 cm

Mathieu-Kocher
TM-11-288 14 cm
TM-11-289 17 cm

TM-11-290 14 cm
TM-11-291 17 cm








TM-11-292 14 cm
TM-11-293 17 cm

Mathieu

Mathieu


Boynton
TM-11-294
12 cm

Gillies with Scissor
TM-11-295 Left handed
TM-11-296 Right handed
16 cm






Needle Holders




Needle Holders with Tungsten Carbide inserts






Micro Needle Holders Stainless Steel



Micro Needle Holders Stainless Steel
with Tungsten Carbide Inserts



Micro Forceps Stainless Steel



Page 15–33

- EXTRACTING FORCEPS



Fig. 1
TM-15-01
Upper Centrals &
Canines / Cuspids



Fig. 2
TM-15-02
Upper Laterals &
Canines / Cuspids



Fig. 3
TM-15-03
Upper Incisors
& Canines



Fig. 4
TM-15-04
Lower Incisors
& Cuspids



Fig. 7
TM-15-07
Upper Premolars / Bicuspid



Fig. 8
TM-15-08
Lower Premolars / Bicuspid



Fig. 13
TM-15-13
Upper Molars, Right



Fig. 17
TM-15-17
Upper Molars, Left



Fig. 18
TM-15-18
Upper Molars



Fig. 18A
TM-15-18A
Upper Molars



Fig. 19
TM-15-19
Upper Third Molars



Fig. 21
TM-15-21
Lower Molars



Extraction Forceps Standard English Pattern



Fig. 22
TM-16-22
Lower Molars



Fig. 23
TM-16-23
Lower Molars
Right



Fig. 24
TM-16-24
Lower Molars
Left




Fig. 29
TM-16-29
Upper Roots




Fig. 29N
TM-16-29N
Upper Roots



Fig. 30
TM-16-30
Upper Roots



Fig. 31
TM-16-31
Upper Roots



Fig. 32
TM-16-32
Lower Molars



Fig. 33
TM-16-33
Lower Roots



Fig. 33A
TM-16-33A
Lower Roots



Fig. 33M
TM-16-33M
Lower Roots



Fig. 41
TM-16-41
Upper Roots




Fig. 44
TM-17-44
Upper Roots



Fig. 45
TM-17-45
Lower Roots



Fig. 46L
TM-17-46L
Lower Roots




Fig. 51
TM-17-51
Upper Roots




Fig. 51A
TM-17-51A
Upper Roots




Fig. 51M
TM-17-51M
Upper Roots



Fig. 52
TM-17-52
Upper Roots



Fig. 52A
TM-17-52A
Upper Teeth
& Resection



Fig. 55
TM-17-55
For Separating
Upper Molars



Fig. 56
TM-17-56
For Separating
Lower Molars



Fig. 67
TM-17-67
Upper Third
Molars




Fig. 67A
TM-17-67A
Upper Third
Molars



Extraction Forceps Standard English Pattern



Fig. 91
TM-18-91
For Separating
Lower Molars



Fig. 99
TM-18-99
Lower Molars



Fig. 99 1/4
TM-18-990
Lower Molars



Fig. 99 1/2
TM-18-991
Lower Molars



Fig. 145
TM-18-145
Very Fine Roots
Lower



Fig. 151 1/2
TM-18-151
Very Fine Roots
Upper



Fig. 349
TM-18-349
Very Fine Roots
Upper



Fig. 359
TM-18-359
Very Fine Roots
Lower



MD 1
TM-18-MD1
Incisors, Bicuspid & Roots
Upper



MD 3
TM-18-MD3
Incisors, Bicuspid & Roots
Lower



MD 2
TM-18-MD2
Upper Molars



MD 4
TM-18-MD4
Lower Molars



Fig. 68
TM-19-68
Lower Roots



Fig. 68A
TM-19-68A
Lower Roots



Fig. 73
TM-19-73
Lower Molars



Fig. 74
TM-19-74
Lower Roots



Fig. 74N
TM-19-74N
Lower Roots



Fig. 75
TM-19-75
Lower Premolars / Bicuspid



Fig. 79
TM-19-79
Lower Third Molars



Fig. 79A
TM-19-79A
Lower Third Molars



Fig. 86A
TM-19-86A
Lower Molars



Fig. 86C
TM-19-86C
Lower Molars



Fig. 87
TM-19-87
Lower Molars



Fig. 89
TM-19-89
Right Upper Molars



Fig. 90
TM-19-90
Left Upper Molars



Extraction Forceps With English Pattern




Fig. 1
TM-20-01
Upper Centrals &
Canines / Cuspids



Fig. 7
TM-20-07
Upper Premolars / Bicuspids



Fig. 13
TM-20-13
Upper Molars, Right



Fig. 17
TM-20-17
Upper Molars, Left



Fig. 18
TM-20-18
Upper Molars



Fig. 22
TM-20-22
Lower Molars



Fig. 33
TM-20-33
Lower Roots



Fig. 51A
TM-20-51A
Upper Roots



Fig. 67A
TM-20-67A
Upper Third
Molars




Fig. 79
TM-20-79
Lower Third Molars



TM-20-10S
Set of 10 pieces



Fig.3
TM-21-03S
Lower Premolars



Fig.5
TM-21-05S
Lower Roots



Fig.6
TM-21-06S
Lower Roots



Fig.7
TM-21-07S
Upper Incisors




Fig.51S
TM-21-51S
Lower Incisors and Canines



Fig.137
TM-21-137S
Upper Molars, left



Fig.139
TM-21-139S
Upper Molars, Right



TM-21-007S
Set of 7 pieces



Extraction Forceps With Anatomical Grip




Fig. S-31
TM-22-S31
Upper Premolars



Fig.34
TM-22-34
Upper Incisors & Cuspids




Fig.35
TM-22-35
Upper Incisors & Cuspids,
Deep Gripping



Fig. 79N
TM-22-79N
Lower Wisdom Teeth




Fig. 46N
TM-22-46N
Lower Premolars



Fig.36
TM-22-36
Lower roots
Incisors & Cuspids



Fig.36N
TM-22-36N
Lower roots
Incisors & Cuspids



Fig.59
TM-22-59
Lower Root Fragments




Fig. 80
TM-23-80
Felsch
Upper Wisdom Teeth
with Retaining Device



Fig. 67L
TM-23-67L
Upper Molars
Left




Fig. 67R
TM-23-67R
Upper Molars
Right




Fig. 67N
TM-23-67N
Upper Molars



Fig. 51B
TM-23-51B
Upper Molars



Fig. 97
TM-23-97
Upper Root Fragments



Fig. 97B
TM-23-97B
Upper Root Fragments



Fig. 122 L
TM-23-122L
Lower Molars & Wisdom Teeth
Left



Fig. 122 R
TM-23-122R
Lower Molars & Wisdom Teeth
Right




Fig. 133 L
TM-23-133L
Lower Very Fine Roots
Left




Fig. 133 R
TM-23-133R
Lower Very Fine Roots
Right



Extraction Forceps "American Pattern"




Standard #1
TM-24-01
For upper incisors & cuspids



Henahan #1A
TM-24-1A
For upper incisors




#2
TM-24-02
Upper & Lower Roots




Woodward #3 FS
TM-24-3FS
For lower molars



#4
TM-24-04
Upper & Lower Premolars




Physick #5
TM-24-05
Lower third molars




Spalitzange #6
TM-24-06
Splitting & Separating Forceps



#7
TM-24-07
Upper Roots



TM-25-08 # 8
Upper & Lower third molars



TM-25-10H # 10 H
For upper Molars, Hook Handle




TM-25-10S # 10S
For upper Molars




TM-25-13 # 13
Upper Incisors




TM-25-14 # 14
Lower Incisors



TM-25-15 # 15
For 1st & 2nd Lower Molars




"Cow Horn"
TM-25-16 # 16
for 1st & 2nd Lower Molars




"Cow Horn"
TM-25-16SK # 16 SK
for children, for Lower Molars




Extraction Forceps "American Pattern"



TM-26-17 # 17
TM-26-17G' # 17 Gold Plated
for 1st & 2nd Lower Molars



TM-26-18L # 18 L
for 1st & 2nd Upper Molars, left




TM-26-18R # 18 R
for 1st & 2nd Upper Molars, right




TM-26-20L # 20 L
Upper Molars left




TM-26-20R # 20 R
Upper Molars right



"Cow Horn"
TM-26-23 # 23
for Lower Molars



TM-26-24 # 24
for Upper Molars




TM-26-27 # 27
for Upper Molars




TM-27-29S # 29 S
Lower teeth for children



TM-27-30S # 30 S
Upper teeth for children



TM-27-32 # 32
for Upper Bicuspids



Parmly
TM-27-32A # 32 A
for upper Bicuspids & Roots



TM-27-34 # 34
Lower Anteriors & Roots



TM-27-35 # 35
Upper Roots



TM-27-36 # 36
Lower Anteriors & Roots




TM-27-37 # 37
Lower Anteriors & Roots




Extraction Forceps "American Pattern"




TM-28-38 # 38
Upper Anteriors & Roots




TM-28-53L # 53 L
TM-28-53LG # 53 L Gold Plated
for upper Molars, left




TM-28-53R # 53 R
TM-28-53RG # 53R Gold Plated
for upper Molars, right




TM-28-62 # 62
for upper & Lower Incisors




TM-28-65 # 65
TM-28-65G # 65 Gold Plated
for upper Incisors, Roots & Fragments




TM-28-69 # 69
TM-28-69G # 69 Gold Plated
for upper & Lower Fragments & Roots




TM-28-85A # 85 A
for Lower Molars & Bicuspids



Nevius
TM-28-88L # 88 L
for Upper Molars, left




Nevius
TM-29-88R # 88 R
For upper Molars, right




TM-29-89 # 89
For upper Molars, right




TM-29-90 # 90
For upper Molars, left




Kells
TM-29-99A # 99 A
Upper Incisors, Canines, Premolars



Kells
TM-29-99C # 99 C
For upper Incisors & Bicuspid



Hull
TM-29-101 # 101
Universal for Bicuspid & Deciduous




TM-29-103 # 103
for Lower Incisors, Bicuspid & Roots




TM-29-104 # 104
Upper third molars




Extraction Forceps "American Pattern"




Nevius
TM-30-110S # 110 S
for Upper Bicuspids & Cuspids




Cryer
TM-30-150 # 150
TM-30-150G # 150 Gold Plated
for Upper Incisors, Bicuspids & Roots




Cryer
TM-30-150A # 150 A
For upper Incisors, Cuspids & Bicuspids



Cryer
TM-30-150AS # 150 AS
For upper Incisors, Cuspids
serrated jaw




Cryer
TM-30-150XAS # 150 XAS
For children, for Upper Incisors &
Cuspids, serrated jaw




Cryer
TM-30-150S # 150 S
TM-30-150SG # 150 S Gold Plated
For children, for Upper Incisors, Bicuspids & Roots




Cryer
TM-30-151 # 151
TM-30-151G # 151 Gold Plated
For Lower Incisors, Bicuspids & Roots




Cryer
TM-30-151A # 151 A
For Lower Incisors, Cuspids
& Bicuspids



Cryer
TM-31-151AS # 151 AS
For Lower Incisors & Cuspids, serrated jaw



TM-31-151XAS # 151 XAS
For children, for Lower Incisors & Cuspids, serrated jaw




Cryer
TM-31-151S # 151 S
TM-31-151SG # 151 S Gold Plated
For children, for Lower incisors, Bicuspids & Roots




TM-31-201 # 201
Upper Incisors, Canines, premolars



TM-31-203 # 203
Lower Incisors & Roots



TM-31-210H # 210 H
For Upper Wisdoms, Hook Handle




TM-31-210S # 210 S
For Upper Wisdoms




TM-31-213 # 213
Upper Centrals & Canines




Extraction Forceps "American Pattern"



TM-32-217 # 217
For Lower Molars




TM-32-222 # 222
For Lower Wisdoms




TM-32-286 # 286
For Upper Incisors, Bicuspid & Roots




TM-32-288 # 288
Lower Canines & Premolars



Cohen
TM-32-310 # 310
Universal Forceps



Nevius
TM-32-320 # 320
Upper Canines



Henahan
TM-32-315 # 315
Upper Incisors




Nevius
TM-32-325 # 325
Upper Canines




Standard
TM-33-330 # 330/4
Upper Third Molars




Woodward
TM-33-335 # 335/FH
Lower Molars




TM-33-17SK
17 SK
For children for Lower Wisdoms



TM-33-23SK
23 SK
For children for Lower Wisdoms



TM-33-150SK
150 SK
For Children, for Upper Incisors,
Bicuspids & Roots



TM-33-151SK
151 SK
For Children, for Lower Incisors,
Bicuspids & Roots



TARIQ MAHMOOD
● ● ● ● ● enterprises

Page 34–73

- MOUTH MIRROR & MIRROR HANDLES
- SCALPEL HANDLES & BLADES
- EXPLORERS SINGLE & DOUBLE ENDED
- PROBES SINGLE ENDED
- GRACEY CURETTES
- PERIODONTAL INSTRUMENTS
- SCALERS
- ENDODONTIC INSTRUMENTS
- FILLING INSTRUMENTS
- CAVITY PREPRATION TIPS & AMALGAM INSTRUMENTS
- AMALGAM PLUGGERS
- BURNISHERS
- COMPOSITE / PLASTIC FILLING INSTRUMENTS
- EXCAVATORS
- CARVERS / FINISHING KNIVES & CLEOID - DISCOID
- WAX MODELLING CARVERS
- CEMENT SPATULAS
- PLASTER KNIVES - ALGINATE & PLASTER SPATULAS
- WAX KNIVES



Fig. 4
TM-34-400
Ø 22 mm



Fig. 5
TM-34-401
Ø 24 mm

- Rhodium coated front surface mirror
- Wide stem offers comfortable cheek retraction
- Non glare satin finish



Fig. 2
18 mm

Fig. 3
20 mm

Fig. 4
22 mm

Fig. 5
24 mm


Pack of 12				
Fig.	Ø mm	Front surface		Rear surface
		Chrome	Rhodium	Magnifying
2	18	-	-	TM-34-406
3	20	-	-	TM-34-410
4	22	TM-34-402	TM-34-404	TM-34-407
5	24	TM-34-403	TM-34-405	TM-34-408



Fig. 000
TM-34-413
Ø 10 mm

Micro mirrors
Rhodium front surface

Fig. 00
TM-34-414
Ø 12 mm



Mirror handles solid



Fig. 1
TM-34-415



Fig. 2
TM-34-416



Fig. 3
TM-34-417



Fig. 4
TM-34-418

Mirror handles hollow



Fig. 1
TM-34-419



Fig. 2
TM-34-420



Fig. 3
TM-34-421



Fig. 4
TM-34-422



Fig. 5
TM-34-423



Fig. 6
TM-34-424



Fig. 7
TM-34-425

Ø 7 mm Hollow Hexagonal 12 cm

Ø 7.5 mm Hollow Hexagonal 14 cm

Ø 7.5 mm Hollow Hexagonal 14 cm

Ø 10 mm Hollow Round 13.5 cm

Ø 9.5 mm Hollow Hexagonal 13.5 cm


Ø 9.5 mm Hollow Hexagonal 13.5 cm

Ø 8mm Soft grip Hollow handle Green




Handles Selection For Scaler

Double Ended



Double Ended

Single Ended





Sterile
Individually Packed
Box of 100 Pieces




Fig. 10 TM-36-444 Fig. 11 TM-36-445 Fig. 12 TM-36-446 Fig. 12D TM-36-447 Fig. 13 TM-36-448 Fig. 15 TM-36-449 Fig. 15C TM-36-450 Fig. 20 TM-36-451 Fig. 22 TM-36-452 Fig. 23 TM-36-453 Fig. 24 TM-36-454 Fig. 25 TM-36-455




Fig. 3 TM-36-456 Fig. 3G TM-36-457 Fig. 4 TM-36-458 Fig. 3L TM-36-459 Fig. 7K TM-36-460 Fig. 7 TM-36-461 TM-36-462 35 mm
TM-36-463 40 mm
TM-36-464 45 mm
TM-36-465 50 mm Fig. 5S TM-36-466 Fig. 5A TM-36-467 Fig. 5 TM-36-468



Micro Scalpel Handles

Micro scalpel & micro mirror handles

Safety Model



TM-37-469 10 cm
TM-37-470 13 cm



TM-37-471
10 cm



TM-37-472
13 cm




TM-37-473



TM-37-474

Scalpel Handle double Blade
15.5 cm

Scalpel Blade removing Forceps



TM-37-475
15 cm



TM-37-476
13.5 cm



TM-37-477
16.5 cm



Scalpel Handle Adjustable
Locks at Angles: 30°, 60°, 90°, 120°, 150°, 180°
TM-37-478



Explorers double ended rigid



Fig. 23/17 TM-38-480



Fig. 23/17A TM-38-481

Explorers double ended rigid



Fig. 23/17 TM-38-482



Fig. 23/17A TM-38-483



Fig. 23/16 TM-38-484

Explorers double ended Flexible



Fig. 2 TM-38-485



Fig. 3 TM-38-486



Fig. 5 TM-38-487



Fig. TU17/23 TM-38-488




Fig. 6/23 TM-38-489




Explorers Double Ended

Expro (explorer-probe combination)




Probes double ended





Explorers single ended rigid

Explorers single ended flexible





Probes Single Ended

Explorers single ended rigid

Explorers single ended rigid



Fig. 2
TM-41-524



Fig. 6
TM-41-525



Fig. 8
TM-41-526



Fig. 9
TM-41-527



Fig. 10
TM-41-528



Fig. 11
TM-41-529



Fig. 12
TM-41-530



Fig. 16
TM-41-531



Fig. 17
TM-41-532



Fig. 23
TM-41-533



Fig. 1
TM-41-534



Fig. 2
TM-41-535



Fig. 3
TM-41-536



Fig. 4
TM-41-537



Fig. 5
TM-41-538



Fig. 6
TM-41-539



Fig. 7
TM-41-540




Fig. 8
TM-41-541



Fig. 9
TM-41-542



Color Coded




TM-42-563
WHO

Endodontic explorers



Solid handle



Endodontic excavators

Solid handle







Fig. 1
TM-44-576

Perforation Repair Instruments




Fig. 1
TM-44-577




Fig. 2
TM-44-578




Fig. 1
TM-44-579




Fig. 2
TM-44-580

Hand Spreaders



Fig. 3
TM-44-581

Solid handle




Fig. 1
TM-44-582




Fig. 2
TM-44-583




Fig. 4
TM-44-584




Fig. 5
TM-44-585




Fig. 6
TM-44-586




Fig. 7
TM-44-587




Fig. 1
TM-44-588




Fig. 2
TM-44-589



Fig. 3
TM-44-590




Fig. 4
TM-44-591



Endodontic Instruments



Fig.2/4
TM-45-592
Molt



Fig.1
TM-45-593
Ruddle
Left



Fig.2
TM-45-594
Ruddle
Right



Fig.1
TM-45-595
Ruddle
Left




Fig.2
TM-45-596
Ruddle
Right


Endodontic Forceps

Micro forceps for removing broken files

Silver points and broken files removing forceps



Endodontic Tweezers






Fig. 1
TM-46-604

Solid handle



Fig. 2
TM-46-605




Fig. 3
TM-46-606




Fig. 4
TM-46-607




Fig. 1
TM-46-608
Luks




Fig. 2
TM-46-609
Luks



Fig. 3
TM-46-610
Luks




Fig. 4
TM-46-611
Luks




Fig. 1
TM-46-612




Fig. 2
TM-46-613




Fig. 3
TM-46-614




Fig. 4
TM-46-615




Fig. 0 TM-47-616




Fig. 1 TM-47-617



Fig. 2 TM-47-618




Fig. 3 TM-47-619



Fig. 4 TM-47-620



Solid handle



Fig. 6 TM-47-621




Fig. 14 TM-47-622



Fig. 17 TM-47-623




Fig. 18 TM-47-624




Fig. 19 TM-47-625




Fig. 17W TM-47-626




Fig. 18W TM-47-627




Fig. 19W TM-47-628



Fig. 0 TM-47-629




Fig. 1 TM-47-630




Fig. 2 TM-47-631




Fig. 3 TM-47-632




Fig. 4 TM-47-633



Solid handle

Palma

Type A



Fig. 1 **TM-48-634**



Fig. 2 **TM-48-635**



Fig. 3 **TM-48-636**



Fig. 4 **TM-48-637**

Type C



Fig. 1 **TM-48-638**



Fig. 2 **TM-48-639**



Fig. 3 **TM-48-640**



Fig. 4 **TM-48-641**



Type B

Palma

Fig. 1 **TM-48-642**




Fig. 2 **TM-48-643**

Fig. 3 **TM-48-644**




Fig. 4 **TM-48-645**



Solid handle

Blacks

Fig. 77/78 **TM-48-646**




Fig. 79/80 **TM-48-647**

Ferrier

Fig. 26 **TM-48-648**




Fig. 27 **TM-48-649**




Fig. 28 **TM-48-650**







Fig. 29 **TM-48-651**




Acron / Wescott



White





Freedman



Solid Handle



Palma






Solid handle














Acron / Wescott



Ladmore




Beavertail / White





Hollow handle


Acron / Wescott



White






Solid handle






Markley




Solid handle




Black





Palma



Palma



Cavity Preparation Tips & Amalgam Instruments



For micro filling



Amalgam well





Pluggers



Fig. 5A TM-54-772



Fig. 3 TM-54-773 Hollenback



Fig. 1 TM-54-774 Super



Fig. OT TM-54-775 Tanner



Fig. 0/1 TM-54-776



Fig. ½ TM-54-777

Burnishers



Fig. 21B TM-54-778 Acom/Wescott



Fig. 1 TM-54-779 Contact Former



Fig. 2 TM-54-780 Contact Former



Fig. 2 TM-54-781 Ladmore



Fig. 3 TM-54-782 Ladmore



Fig. 27/29 TM-54-783



Fig. 1 TM-54-784



Fig. 1L TM-54-785



Fig. 155 TM-54-786



Fig. 156(A6) TM-55-787



Fig. 8A TM-55-788



Fig. 1 TM-55-789 Boghosian



Fig. 2 TM-55-790 Boghosian



Fig. 1 TM-55-791



Fig. 2 TM-55-792



Fig. 1 TM-55-793 Goldstein



Fig. 2 TM-55-794 Goldstein



Fig. 3 TM-55-795 Goldstein



Fig. 4 TM-55-796 Goldstein



Fig. 1 TM-55-797 Goldstein-Mini



Fig. 3 TM-55-798 Goldstein-Mini



Fig. 4 TM-55-799 Goldstein-Mini



Fig. 1 TM-55-800 Goldstein-Micro



Fig. 4/5 TM-55-801 Gregg



Fig. 0 TM-55-802 Heidemann



Fig. 1 TM-55-803 Heidemann



Fig. 2 TM-55-804 Heidemann



Fig. W3 TM-55-805



Fig. 1 TM-55-806



Fig. 6 TM-55-807



Fig. 178 TM-55-808



Fig. 49 TM-55-809



Fig. 179 TM-55-810



Fig. 6 TM-56-811




Fig. 178 TM-56-812




Fig. 179 TM-56-813




Fig. 5/7 TM-56-814

Gregg




Fig. 4/5 TM-56-815

Woodson




Fig. 1 TM-56-816



Fig. 2 TM-56-817




Fig. 3 TM-56-818

Heidemann



Fig. 0 TM-56-819




Fig. 1 TM-56-820




Fig. 2 TM-56-821




Fig. 0 TM-56-822



Fig. 1 TM-56-823




Fig. 2 TM-56-824



Fig. 3 TM-56-825




Fig. 4 TM-56-826




Fig. 5 TM-56-827




Fig. 6 TM-56-828




Fig. 7 TM-56-829




Fig. 8 TM-56-830




Fig. 9 TM-56-831




Fig. 10 TM-56-832




Fig. 11 TM-56-833



























Fig. 12 TM-56-834

Plastic Filling Instruments

Hollow handle



Solid handle






Fig. 1 TM-58-860




Fig. 6 TM-58-861



Fig. 178 TM-58-862



Fig. 179 TM-58-863

Heidemann




Fig. 0 TM-58-864




Fig. 1 TM-58-865




Fig. 2 TM-58-866




Fig. 7 TM-58-867



Fig. 8 TM-58-868



Fig. 9 TM-58-869



Fig. 10 TM-58-870




Fig. 11 TM-58-871




Fig. 12 TM-58-872

Gregg




Fig. 4/5 TM-58-873


Applicators for calcium hydroxide



TM-58-874



TM-58-875



TM-58-876




TM-58-877

Applicators for medicaments



TM-58-878
Orban





Solid handle




Fig. 204
TM-60-892




Fig. CI 2/3
TM-60-899




Fig. 204S
TM-60-893




Fig. 6
TM-60-900
Crane-Kaplan



Fig. 204SD
TM-60-894




Fig. 6/7
TM-60-901
Hygienist




Fig. 31/32
TM-60-895
Jacquette




Fig. UI5/33
TM-60-902
Jacquette/Towner




Fig. 30/33
TM-60-896
Jacquette




Fig. 33/H5
TM-60-897
Jacquette




Fig. 2/3
TM-60-903
Taylor




Fig. UI5/30
TM-60-898
Jacquette/Towner




Fig. 7/8
TM-60-904
Younger-Good




Fig. 204
TM-61-905




Fig. CI 2/3
TM-61-912




Fig. 204S
TM-61-906




Fig. 6
TM-61-913
Crane-Kaplan



Fig. 204SD
TM-61-907




Fig. 6/7
TM-61-914
Hygienist




Fig. 31/32
TM-61-908
Jacquette



Fig. U15/33
TM-61-915
Jacquette/Towner




Fig. 30/33
TM-61-909
Jacquette




Fig. 33/H5
TM-61-910
Jacquette




Fig. 2/3
TM-61-916
Taylor



Fig. U15/30
TM-61-911
Jacquette/Towner




Fig. 7/8
TM-61-917
Younger-Good




Solid Handle



Darby Perry



Standard













Standard










Palma



Hygienist




Jacquette




Hollow Handle




Hygienist











Jacquette





Standard



Standard



Towner



Periodontal chisels



Fig. 1 TM-63-962 TG



Fig. 2 TM-63-963 TG

Universal Curettes Columbia

Solid handle



Hollow handle



Solid handle

Hollow handle

Universal

Fig. 13R/14L

TM-63-964

TM-63-967




Anterior teeth

Fig. 2R/2L

TM-63-965

TM-63-968



Posterior teeth

Fig. 4R/4L

TM-63-966

TM-63-969

Hoe scalars

Solid handle



Hollow handle



Goldman-Fox



Solid handle

Hollow handle

Fig. 5

TM-63-970

TM-63-972




Fig. 6

TM-63-971

TM-63-973



Solid handle



Hollow handle



Anterior teeth




Fig. 1/2

Anterior teeth




Fig. 3/4

Anterior and
Premolars




Fig. 5/6

Premolars and
Molars




Fig. 7/8

Molars buccal and
lingual surfaces



Fig. 9/10

Posterior mesial
surfaces




Fig. 11/12

Posterior distal
surfaces




Fig. 13/14

Solid handle

Fig. 1/2	TM-64-974
Fig. 3/4	TM-64-975
Fig. 5/6	TM-64-976
Fig. 7/8	TM-64-977
Fig. 9/10	TM-64-978
Fig. 11/12	TM-64-979
Fig. 13/14	TM-64-980

TM-64-980-S

Gracey set of 7 pieces

Hollow handle

Fig. 1/2	TM-64-981
Fig. 3/4	TM-64-982
Fig. 5/6	TM-64-983
Fig. 7/8	TM-64-984
Fig. 9/10	TM-64-985
Fig. 11/12	TM-64-986
Fig. 13/14	TM-64-987

TM-64-987-S

Gracey set of 7 pieces



Fig. 1/2
TM-65-988
Grey
Anterior teeth



Fig. 3/4
TM-65-989
Black
Anterior teeth



Fig. 5/6
TM-65-990
Yellow
Anterior and Premolars



Fig. 7/8
TM-65-991
Green
Premolars and Molars



Fig. 9/10
TM-65-992
Red
Molars buccal and lingual surfaces



Fig. 11/12
TM-65-993
Orange
Posterior mesial surfaces



Fig. 13/14
TM-65-994
Blue
Posterior distal surfaces

TM-65-994-S
Set of 7 pieces
Fig. 1/2, 3/4, 5/6, 7/8, 9/10, 11/12, 13/14



Fig. 15/16
TM-65-995
Coffee
Posterior mesial surfaces




Fig. 17/18
TM-65-996
Light Blue
Posterior distal surfaces




Fig. 11/14
TM-65-997
Blue
Mesial/Distal surfaces



Fig. 12/13
TM-65-998
Orange
Mesial/Distal surfaces



- Ergonomic, large diameter silicone handle for better grip
- Ultra light weight
- Special grade stainless steel tips ensure longer life
- Rigid pattern, thicker tips than standard Gracey



TM-66-1006-S
Set of 7 pieces / Rigid

TM-66-1013-S
Set of 7 pieces / Standard




Fig.3 TM-67-1014 Hollenback




Fig.2/3 TM-67-1018 Frahm



Fig.1 TM-67-1015 Interproximal



Fig.1 TM-67-1019 Wards




Fig.1 TM-67-1016 Interproximal



Fig.2 TM-67-1020 Wards




Fig.2 TM-67-1017 Interproximal




Fig.8 TM-67-1021 Wiland



Hollow handle



Fig.6T TM-68-1022 Tanner



Fig.2 TM-68-1023 USC



Fig.2S TM-68-1024 USC

Cleoid - Discoid



Solid handle




Fig.1 TM-68-1025




Fig.3 TM-68-1027




Fig.2 TM-68-1026




Fig.4 TM-68-1028

Cleoid - Discoid




Fig.4/5 TM-68-1029




Fig.89/92 TM-68-1031




Fig.3/6 TM-68-1030




Fig.90B TM-68-1032

Hollow Handle



Fig.1 TM-69-1033



Fig.2 TM-69-1034



Fig.4 TM-69-1036



Fig.3 TM-69-1035



Fig.5 TM-69-1037

Solid Handle



Fig.1 TM-69-1038



Fig.2 TM-69-1039



TM-69-1043-S
P.k. Thomas
Set of 5 Pieces



Fig.4 TM-69-1041



Fig.3 TM-69-1040



Fig.5 TM-69-1042

Hollow Handle



Fig.7 TM-69-1044 Beale



Fig.31 TM-69-1045 Gritman



Fig.10 TM-69-1047 Le cron



Fig.10 TM-69-1046 Hylin



Fig.10 TM-69-1048 Zahle

Solid Handle



Fig.1 TM-69-1049 Vehe



Fig.31 TM-69-1050 Gritman



Fig.10 TM-69-1052



Fig.7 TM-69-1053 Beale



Fig.12 TM-69-1051



Fig.1 TM-69-1054 Roach



Fig. 1 13 cm



Fig. 2 14 cm



Fig. 3 14 cm



Fig. 4 13.5 cm



Fig. 5 14 cm



Fig. 6 11.7 cm



Fig. 7 15 cm



Fig. 1 Zahle 12 cm



Fig. 2 Zahle 12.5 cm



Fig. 3 Zahle 12.5 cm



Fig. 4 Zahle 12 cm



Fig. 5 Zahle 12 cm



Fig. 6 Zahle 11.5 cm

TM-70-1055

TM-70-1056

TM-70-1057

TM-70-1058

TM-70-1059

TM-70-1060

TM-70-1061

TM-70-1062

TM-70-1063

TM-70-1064

TM-70-1065

TM-70-1066

TM-70-1067



Fig. 1 Hylin 13.5 cm



Fig. 3 Hylin-Ceramica 13.5 cm



Fig. 5 Hylin 13.5 cm



Fig. 6 Hylin-Lecron 13.5 cm



Fig. 1 Lecron-mini 10 cm



Fig. 2 Lecron-mini 13 cm



Fig. 3 Lecron-mini 13 cm



Fig. 1 Lecron 16.5 cm



Fig. 2 Lecron 16.5 cm



Fig. 3 Lecron 16.5 cm



Fig. 4 Lecron 16.5 cm



Fig. 7 Lecron 16 cm



Fig. 8 Lecron 16 cm

TM-70-1068

TM-70-1069

TM-70-1070

TM-70-1071

TM-70-1072

TM-70-1073

TM-70-1074

TM-70-1075

TM-70-1076

TM-70-1077

TM-70-1078

TM-70-1079

TM-70-1080



Fig.1 TM-71-1081



Fig.2 TM-71-1082



Fig.3 TM-71-1083



Fig.1 TM-71-1084



Fig.2 TM-71-1085



Fig.3 TM-71-1086



Fig.00 TM-71-1087



Fig.0 TM-71-1088




Fig.1 TM-71-1089




Fig.2 TM-71-1090



Fig.3 TM-71-1091



Fig.4 TM-71-1092



Weston TM-71-1093



TM-71-1094



Plaster knives with wooden handle



TM-72-1106 Fahnenstock 13 cm



TM-72-1108 Lessmann 13 cm



TM-72-1107 Fahnenstock 17.5 cm



TM-72-1109 Lessmann 17.5 cm

Wax knives with wooden handle



TM-72-1110 Fahnenstock 13 cm



TM-72-1112 Lessmann 13 cm



TM-72-1111 Fahnenstock 17.5 cm



TM-72-1113 Lessmann 17.5 cm



TM-72-1114 Con Scodelino 13 cm



TM-72-1115 Con Scodelino 17.5 cm

Wax knives Ermert



TM-72-1116 Wooden Handle



TM-72-1117 Plastic Handle



TM-72-1118 Metal Handle

Wax knives with metal handle



TM-72-1119 Fehnenstock 13 cm



TM-72-1120 Fehnenstock 17.5 cm



TM-72-1121 Lessmann 13 cm




TM-72-1122 Lessmann 17.5 cm

Plaster Spatulas

Flexible Spatulas



Rigid Spatulas





TARIQ MAHMOOD
● ● ● ● ● enterprises

Page 74–83

o ROOT ELEVATORS



Root Elevators

Luxators



Fig. 1
TM-74-1141



Fig. 2
TM-74-1142



Fig. 3
TM-74-1143



Fig. 4
TM-74-1144



Fig. 5
TM-74-1145



Fig. 6
TM-74-1146



Fig. 7
TM-74-1147



TM-74-1150-S
Set of 7 Pcs

Bein



Fig. 1
TM-74-1151



Fig. 2
TM-74-1152



Fig. 3
TM-74-1153



Fig. 4
TM-74-1154



Fig. 5
TM-74-1155



Fig. 6
TM-74-1156

Apical



Fig. 301
TM-74-1157



Fig. 302
TM-74-1158
Left




Fig. 303
TM-74-1159
Right



Fig. 304
TM-74-1160



Seldin



Fig. 34
TM-75-1161



Fig. 34S
TM-75-1162



Fig. 1L
TM-75-1163
Left



Fig. 1R
TM-75-1164
Right



Fig. 4L
TM-75-1165
Left



Fig. 4R
TM-75-1166
Right

Cryer



Fig. 25
TM-75-1167
Left



Fig. 26
TM-75-1168
Right



Fig. 31
TM-75-1169
Left



Fig. 32
TM-75-1170
Right



Fig. 44
TM-75-1171
Left



Fig. 45
TM-75-1172
Right

Coupland




Fig. 1
TM-75-1173




Fig. 2
TM-75-1174




Fig. 3
TM-75-1175



Fig. 4
TM-75-1176
Left



Fig. 5
TM-75-1177
Right

Schmeckebier





Root Elevators Bien



Fig. 1
TM-76-1178



Fig. 2
TM-76-1179



Fig. 3
TM-76-1180



Fig. 4
TM-76-1181



Fig. 5
TM-76-1182



Fig. 6
TM-76-1183



Fig. 7
TM-76-1184



Fig. 8
TM-76-1185



Fig. 9
TM-76-1186



Fig. 10
TM-76-1187



Fig. 1
TM-76-1188



Fig. 2
TM-76-1189



Fig. 3
TM-76-1190



Fig. 1
TM-76-1191



Fig. 2
TM-76-1192



Fig. 3
TM-76-1193



Fig. 1
TM-76-1194
Bayonet Distal



Fig. 2
TM-76-1195
Bayonet Distal



Fig. 3
TM-76-1196
Bayonet Distal



Fig. 1
TM-76-1197
Bayonet Mesial



Fig. 2
TM-76-1198
Bayonet Mesial



Fig. 3
TM-76-1199
Bayonet Mesial



Flohr



Fig. 1
TM-77-1200



Fig. 2
TM-77-1201
Left



Fig. 3
TM-77-1202
Right

Friedman



Fig. 1
TM-77-1203




Fig. 2
TM-77-1204
Left



Fig. 3
TM-77-1205
Right

Lindo-Levian



Fig. LLL
TM-77-1206



Fig. LLM
TM-77-1207




Fig. LLS
TM-77-1208




Fig. LLL
TM-77-1209




Fig. LLM
TM-77-1210




Fig. LLS
TM-77-1211



Fig. LLM
TM-77-1212
Bayonet Distal




Fig. LLM
TM-77-1213
Bayonet Mesial

Hylin



Fig. 1
TM-77-1214



Fig. 2
TM-77-1215



Fig. 3
TM-77-1216
Left




Fig. 4
TM-77-1217
Right

Bein



Fig. 1
TM-77-1218



Fig. 6
TM-77-1219




Fig. 7
TM-77-1220

Flohr



Fig. 1
TM-77-1221



Fig. 2
TM-77-1222
Left




Fig. 3
TM-77-1223
Right



Root Elevators

Schlemmer



Fig. 1
TM-78-1224
Left



Fig. 2
TM-78-1225
Right

Coleman / London Hospital



Fig. 20
TM-78-1228
Coleman



Fig. 21
TM-78-1229
London Hospital
Left



Fig. 22
TM-78-1230
London Hospital
Right

Schlemmer



Fig. 1
TM-78-1226
Left



Fig. 2
TM-78-1227
Right

Warwick-James



TM-78-1231
Straight



TM-78-1232
Left



TM-78-1233
Right

Root Tip Picks / Heidebrink



Fig. 2
TM-78-1235
Left



Fig. 1
TM-78-1234



Fig. 3
TM-78-1236
Right



Fig. 13/14
TM-78-1237
Double ended



Fig. 1
TM-78-1238



Fig. 2
TM-78-1239
Left




Fig. 3
TM-78-1240
Right



TM-79-1241
Lecluse Fig. 1



Potts



Fig. 6
TM-79-1241
Left



Fig. 7
TM-79-1242
Right

Winter



Fig. 1L
TM-79-1244
Left



Fig. 1R
TM-79-1245
Right



Fig. 11L
TM-79-1246
Left



Fig. 11R
TM-79-1247
Right



Fig. 12L
TM-79-1248
Left



Fig. 12R
TM-79-1249
Right



Fig. 13L
TM-79-1250
Left



Fig. 13R
TM-79-1251
Right




Fig. 13L
TM-79-1252
Left



Fig. 13R
TM-79-1253
Right

Walter F. Barry



Fig. 320
TM-79-1254
Left



Fig. 321
TM-79-1255
Right



Fig. 322
TM-79-1256
Left



Fig. 323
TM-79-1257
Right



Fig. 350
TM-79-1258
Left



Fig. 351
TM-79-1259
Right



Root Elevators



Fig. 1
TM-80-1260



Fig. 2
TM-80-1261



Fig. 3
TM-80-1262



Fig. 44
TM-80-1263
Left



Fig. 45
TM-80-1264
Right



Fig. 2
TM-80-1265
Left



Fig. 3
TM-80-1266
Right



TM-80-1267-S
Set of 7 Pcs



Fig. 1
TM-80-1268



Fig. 2
TM-80-1269



Fig. 3
TM-80-1270



Fig. 44
TM-80-1271
Left



Fig. 45
TM-80-1272
Right



Fig. 2
TM-80-1273
Left



Fig. 3
TM-80-1274
Right



TM-80-1275-S
Set of 7 Pcs



Fig. 1
TM-81-1276



Fig. 2
TM-81-1277



Fig. 3
TM-81-1278



Fig. 44
TM-81-1279
Left



Fig. 45
TM-81-1280
Right



Fig. 2
TM-81-1281
Left



Fig. 3
TM-81-1282
Right



TM-81-1283-S
Set of 7 Pcs

Chompret



Fig. 1
TM-81-1284



Fig. 2
TM-81-1285



Fig. 3
TM-81-1286



Fig. 4
TM-81-1287





Root Elevators

Serrated



Fig. 1
TM-82-1288



Fig. 2
TM-82-1289



Fig. 1
TM-82-1290



Fig. 3
TM-82-1291



Fig. 4
TM-82-1292



Fig. 6
TM-82-1293



Fig. 7
TM-82-1294

Serrated Elevators TM-82-1295-S Set of 7 Pcs

Luxators



Fig. 1
TM-82-1296



Fig. 2
TM-82-1297



Fig. 3
TM-82-1298



Fig. 4
TM-82-1299



Fig. 5
TM-82-1300



Fig. 6
TM-82-1301




Fig. 7
TM-82-1302

Luxators TM-82-1303-S Set of 7 Pcs






Fig. 2
TM-83-1304
Flohr



Fig. 3
TM-83-1305
Flohr



Fig. 1
TM-83-1306
Bein




Fig. 6
TM-83-1307
Bein



Fig. 7
TM-83-1308
Bein

TM-83-1309-S Set of 5 Pcs

British Models




Fig. 1
TM-83-1310
Coupland



Fig. 2
TM-83-1311
Coupland



Fig. 3
TM-83-1312
Coupland



TM-83-1313
Warwick-James



TM-83-1314
Warwick-James



TM-83-1315
Warwick-James



Fig. 30
TM-83-1316
Cryer



Fig. 31
TM-83-1317
Cryer

TM-83-1318-S Set of 8 Pcs

Root Tip Elevators



Fig. 1
TM-83-1319
Heidebrink




Fig. 2
TM-83-1320
Heidebrink




Fig. 3
TM-83-1321
Heidebrink



Fig. 11
TM-83-1322
Bernhard



Fig. 6
TM-83-1323
Bernhard



Fig. 7
TM-83-1324
Bernhard

TM-83-1324-S Set of 6 Pcs



TARIQ MAHMOOD
enterprises

Page 84–110


- MOUTH GAGS
- TONGUE DEPRESSORS
- LIPS & CHEEK RETRACTORS
- RETRACTORS & FORCEPS
- BONE INSTRUMENTS
- IMPLANTOLOGY INSTRUMENTS
- PERIOSTEAL ELEVATORS
- PERIODONTAL INSTRUMENTS
- BONE CURETTES & FILES
- SURGERY INSTRUMENTS



Molt
TM-84-1325 10 cm
TM-84-1326 11 cm
TM-84-1327 14 cm



Roser-König
TM-84-1328 16 cm
TM-84-1329 19 cm



TM-84-1330
Denhardt
13 cm



TM-84-1331
Fergusson - Ackland
14.5 cm




TM-84-1332
Heister
Maunders Wood
9.5 cm



TM-84-1333
Obwegeser 11 cm
Spreading 47 mm



TM-84-1334
Heister
13 cm



Davis-Boyle Tongue Depressors
With fixed suction tube

Fig. 1
TM-85-1335
110 x 28 mm



Fig. 2
TM-85-1336
95 x 27 mm



Fig. 3
TM-85-1337
80 x 27 mm



Fig. 4
TM-85-1338
65 x 26 mm



Fig. 5
TM-85-1339
30 x 20 mm



TM-85-1340
Davis-Boyle
Complete, consisting of frame
& 5 tongue depressors No. 1 to 5



TM-86-1341
Mc-Ivor
Frame for Mouth Gag



Mc-Ivor
Tongue Depressors

Fig. 1
TM-86-1342
100 x 25 mm

Fig. 2
TM-86-1343
85 x 25 mm

Fig. 3
TM-86-1344
75 x 25 mm

TM-86-1345
Mc-Ivor
Complete, with 3 tongue depressors



TM-86-1346
Dingmann
Complete, with frame, 3 tongue depressors,
& 2 adjustable lateral blades



Dingmann
Tongue Depressors
13.5 cm

Fig. 1
TM-86-1347
62 x 26 mm

Fig. 2
TM-86-1348
65 x 30 mm

Fig. 3
TM-86-1349
78 x 32 mm



Jennings

TM-87-1350 10.5 cm
TM-87-1351 12.5 cm
TM-87-1352 14.5 cm
TM-87-1353 16.5 cm




Whitehead Mouth Gag

TM-87-1354 11 cm
TM-87-1355 13 cm
TM-87-1356 15 cm



Cheek Retractor

TM-87-1357 Small
TM-87-1358 Medium
TM-87-1359 Large



TM-88-1360
Hartmann
15 cm



TM-88-1361
Doyen Tongue Spatula
13.5 cm




TM-88-1362
14.5 cm




TM-88-1363
18 cm




TM-88-1364
Bruenings
19 cm




TM-88-1365
Tobold
15 cm




TM-88-1366
Tobold
16 cm



TM-88-1367
14 cm



TM-88-1368
15 cm





Lips & Cheek Retractors



TM-89-1369
Obwegeser Chin
Retractor
16.5 cm



Obwegeser 16 cm
TM-89-1370 8 mm
TM-89-1371 10 mm



TM-89-1372 43 x 11 mm
TM-89-1373 57 x 11 mm
TM-89-1374 83 x 16 mm



Obwegeser
TM-89-1375
Spina-nasal's retractor
22 cm



Obwegeser
TM-89-1376
Ramus retractor for
separation of ramus & tissue
22 cm



TM-89-1377
Simplex
10.5 cm



TM-89-1378
Kilner
15 cm



TM-89-1379
Minnesota
14 cm



TM-89-1380
Cawood-Minnesota
15 cm



TM-89-1381
Austin



TM-90-1382
Hilger
Fits any simple stem mouth mirror handle



TM-90-1383
Hilger
with handle



TM-90-1384
20 cm

Finger protectors



TM-90-1385 Small
TM-90-1386 Medium
TM-90-1387 Large

Langenbeck



TM-90-1388
Oringer
11.5 cm



TM-90-1389
Obwegeser Cheek Retractor
12cm



TM-90-1390
Columbia
14.5 cm



Sternberg


TM-90-1391 12 cm
TM-90-1392 14 cm
TM-90-1393 15 cm



Vestibulum
TM-90-1394




TM-90-1395




TM-90-1396

Implants Retractors

Retractors & Forceps



Senn-Miller
16 cm



16 cm



16 cm



TM-91-1407
Collin Tongue holding forcep
17 cm



Young Tongue holding forcep
TM-91-1408 15 cm
TM-91-1409 17 cm




TM-91-1411
Farabeuf 12 cm
Pair



TM-91-1412
Parker-Langenbeck 21 cm
Pair




Kocher
21.5 cm



TM-92-1413
55 x 11 mm

TM-92-1414
40 x 11 mm



TM-92-1415
55 x 11 mm

TM-92-1416
40 x 11 mm




Middeldorph



TM-92-1417
14 x 17 mm
21.5 cm

TM-92-1418
20 x 22 mm
22 cm




TM-92-1419
26 x 30 mm
23.5 cm




■ Langenbeck
22 cm




● Langenbeck
22 cm




■ TM-92-1420
● TM-92-1421




■ TM-92-1422
● TM-92-1423



■ TM-92-1424
● TM-92-1425



■ TM-92-1426
● TM-92-1427





■ TM-92-1428
● TM-92-1429




Bone Instruments


Chisels & Gouges




Osteotomes



Pocket Markers



Elevators & Raspatories





PDL Elevators



Instruments
Tray



TM-94-1457 Distal, 2.5 mm Micro Edge



TM-94-1458 Mesial, 2.5 mm Micro Edge



TM-94-1459 Distal, 4 mm Micro Edge



TM-94-1460 Mesial 4 mm Micro Edge



TM-94-1461 2.5 mm Micro Edge, Straight



TM-94-1462 2.5 mm Micro Edge, Curved



TM-94-1463 3.4 mm Straight

TM-94-1463-S Set Of 7 Pieces



TM-94-1464 Distal, 2.5 mm Micro Edge



TM-94-1465 Distal, 2.5 mm Micro Edge



TM-94-1466 Mesial, 2.5 mm Micro Edge



TM-94-1467 Mesial 4 mm Micro Edge



TM-94-1468 2.5 mm Micro Edge, Straight





TM-94-1469 2.5 mm Micro Edge, Curved



TM-94-1470 3.4 mm Straight

TM-94-1470-S Set Of 7 Pieces

Osteotomes





TM-96-1485
Bone Collector Chisel Straight



TM-96-1487
Straight blade for bone collector



TM-96-1486
Bone Collector Chisel Curved



TM-96-1488
Curved blade for bone collector



TM-96-1489
Ratchet Screw driver handle
Hexagonal head



TM-96-1490
Bone Morselizer with Titanium Grinding



TM-96-1491
Bone Mixing Basin
Dia 35 mm




TM-96-1492
Manual Easy Bone Mill



TM-96-1493
Bone Craker Titanium



Sinus Lift Currettes





Sinus Lift Elevators

Hollow handle



Hollow handle



TiN Coated



TM-98-1510
Fig. 1



TM-98-1518
Fig. 1



TM-98-1511
Fig. 2



TM-98-1519
Fig. 2



TM-98-1512
Fig. 3



TM-98-1520
Fig. 3



TM-98-1513
Fig. 4



TM-98-1521
Fig. 4



TM-98-1514
Fig. 5



TM-98-1522
Fig. 5



TM-98-1515
Fig. 6



TM-98-1523
Heights Fig. 6



TM-98-1516
Fig. 7



TM-98-1524
Fig. 7



TM-98-1517
Fig. 8



TM-98-1525
Fig. 8


Bone Collector



TM-99-1526

TM-99-1527
Disposable stainless steel cartridge
Pack of 10 pieces

Applicators for bone grafting material



TIN coated

TM-99-1528 Applicator



TM-99-1529 Applicator & compactor

Compactors



TM-99-1530 Hollow Handle



TM-99-1531 Hollow Handle

TIN coated



TM-99-1532 Hollow Handle



TM-99-1533 Hollow Handle

Membrane placement instruments




TM-99-1534




TM-99-1535
11 cm



Osteotomes



Dilator



Spreader

A-Form
Ø 2 mm Ø 3 mm Ø 3.3 mm Ø 4 mm



14.5 mm



TM-100-1536
14.5 cm



TM-100-1537
14.5 cm



TM-100-1538
14.5 cm



TM-100-1539
14.5 cm

D-Form
Convex Tips
Ø 3 mm Ø 4 mm Ø 5 mm



TM-100-1540
14.5 cm




TM-100-1541
14.5 cm



TM-100-1542
14.5 cm



FREER Osteotomos



FREER Osteotomos

Ø 1.6 mm Ø 1.9 mm Ø 2.8 mm Ø 3.1 mm Ø 3.9 mm



TM-100-1543



TM-100-1544



TM-100-1545



TM-100-1546



TM-100-1547

Ø 3.0 mm Ø 3.3 mm Ø 4.0 mm



TM-100-1548




TM-100-1549



TM-100-1550



Osteotomes





TM-102-1565
Interior 2 mm
Exterior 3 mm



TM-102-1566
Interior 3 mm
Exterior 4 mm



TM-102-1567
Interior 4 mm
Exterior 5 mm



TM-102-1568
Interior 5 mm
Exterior 6 mm



TM-102-1569
Interior 6 mm
Exterior 7 mm



TM-102-1570
Interior 7 mm
Exterior 8 mm



TM-102-1571
Interior 8 mm
Exterior 9 mm



TM-102-1572
Interior 9 mm
Exterior 10 mm



TM-102-1573
Trephine Drills Stand
Complete with drills



TM-102-1574
Trephine Drills Stand
Complete with drills

Tissue Punches



Keyes



TM-102-1575
3 mm



TM-102-1576
4 mm



TM-102-1577
5 mm



TM-102-1578
8 mm



Keyes



TM-102-1579
3 mm



TM-102-1580
4 mm



TM-102-1581
5 mm



TM-102-1582
6 mm

Bone Thickness Gauges



Fig. 1
Castroviejo
TM-103-1583



Fig. 2
Castroviejo
TM-103-1584



Fig. 3
Castroviejo
TM-103-1585



Fig. 4
Castroviejo
TM-103-1586



Fig. 5
Castroviejo
TM-103-1587



TM-103-1588
Bone Thickness Gauge



TM-103-1589
Bone Caliper
15 cm



TM-103-1590
Sinus Rongeur
Delicate bone cutting forceps
for opening the sinus window
4 mm bite, TiN coated

Bone Injectors



Fig. 1
TM-103-1591



Fig. 2
TM-103-1592



Fig. 3
TM-103-1593




Fig. 4
TM-103-1594



TiN coated bone block clamps

To hold & stabilize block graft against recipient site while placing screws.
Slotted tips act as a guide for precise drilling.



TM-103-1595
Posterior

TM-103-1596
Anterior

Depth measuring probes



TM-103-1597
Titanium tip, silicone handle & laser marked



Periotome



Fig. 1 **TM-104-1598**

Implant periosteal elevators & raspatories



Fig.1 **TM-104-1600**



Fig.2 **TM-104-1601 Sharp-Blunt**



Fig.1 **TM-104-1602** Atraumatic papilla elevator

Mallets



TM-104-1603
Partsch
Stainless Steel head
19 cm



TM-104-1604
Partsch
Lead filled head
19 cm



TM-104-1605
Mead 19 cm
TM-104-1606
Spare plastic disks, pair



TM-104-1607
18.5 cm
TM-104-1608
Spare plastic disks, pair



TM-104-1609
19 cm
TM-104-1610
Spare plastic disks, pair



TM-104-1611
20 cm
TM-104-1612
Spare plastic disks, pair



Fig.15 TM-105-1613 Freer



Fig.16 TM-105-1614 Freer




Fig.14 TM-105-1615 FGoldman-Fox




Fig.20 TM-105-1616 Hirschfield



Fig.1 TM-105-1617 Hopkins



Fig.2 TM-105-1618 Hopkins



Fig.3 TM-105-1619 Mead/West



Fig.9 TM-105-1620 Molt



Fig.3 TM-105-1621 Prichard



Fig.9 TM-105-1622 Molt



Fig.1 TM-106-1623 Buser Blue



Fig.3 TM-106-1624 Prichard Blue

Implant periosteal elevators



Fig.1 TM-106-1625 Buser



Fig.9A TM-106-1626 Allen



Fig.25 TM-106-1627 17 cm

Gingivectomy Knives

Hollow handle




Fig.15/16 TM-106-1628 Kirkland



Fig.1/2 TM-106-1629 Orban

Solid handle



Fig.15/16 TM-106-1630 Kirkland




Fig.1/2 TM-106-1631 Orban











Fig.1 TM-106-1632 Orban





Fig.2 TM-106-1633 Orban





Bone Curettes





	Solid handle	Hollow handle
Fig. 1	TM-108-1647	TM-108-1651
Fig. 2	TM-108-1648	TM-108-1652
Fig. 3	TM-108-1649	TM-108-1653
Fig. 4	TM-108-1650	TM-108-1654



	Solid handle	Hollow handle
Fig. 85	TM-108-1655	TM-108-1659
Fig. 86	TM-108-1656	TM-108-1660
Fig. 87	TM-108-1657	TM-108-1661
Fig. 88	TM-108-1658	TM-108-1662






Bone Currettes & Files




Bone currettes



Molt
16.5 cm



Williger
14 cm



Bone Files



Fig. 10
TM-109-1671



Fig. 12A
TM-109-1672



Fig. 33
TM-109-1673



Fig. 12
TM-109-1674
Howard



Fig. 45
TM-109-1675
Miller




Fig. 64
TM-109-1676
Miller







Skin Hooks



Probes




Self-Retaining Retractors



Suture Forceps



Suction Tubes





TARIQ MAHMOOD
enterprises

Page 111–115

- SELF & MANUAL ASPIRATING SYRINGES
- INTRALIGAMENTAL SYRINGES



TM-111-1700
Type - A
Metric Thread
1.8 ml



TM-111-1701
Type - C
Metric Thread
1.8 ml


Manual Aspirating Syringes



TM-111-1702
Type - A
Metric Thread
1.8 ml



TM-111-1703
Type - C
Metric Thread
1.8 ml



TM-111-1704
P - Ject
Metric Thread
1.8 ml



TM-111-1705
Type - A
Metric Thread
1.8 ml



TM-111-1706
Euro
Metric Thread
1.8 ml




Manual Aspirating Syringes



TM-112-1707
Euro
Metric Thread
1.8 ml




TM-112-1708
T - Handle
Metric Thread
1.8 ml



TM-112-1709
3 - Ring
Metric Thread
1.8 ml



TM-112-1710
Euro breech loading, Ring handle
Metric Thread
1.8 ml



TM-112-1711
Euro breech loading, T - handle
Metric Thread
1.8 ml



TM-113-1712
Euro
Metric Thread
1.8 ml




TM-113-1713
3 - Ring
Metric Thread
1.8 ml



TM-113-1714
2 - Ring
Metric Thread
1.8 ml



TM-113-1715
Euro breech loading, Ring handle
Metric Thread
1.8 ml



TM-113-1716
Euro breech loading, T - handle
Metric Thread
1.8 ml



TM-114-1717

Curved
Metric Thread
1.8 ml



TM-114-1718

Curved
Metric Thread
1.8 ml



TM-114-1719

Straight
Metric Thread
1.8 ml



TM-114-1720

Swedish Type
Metric Thread
1.8 ml





TARIQ MAHMOOD
enterprises

Page 116–129

- RUBBER DAM INSTRUMENTS
- MATRIX RETAINERS & RUBBER DAM FRAMES
- LIGATURE & SEPARATORS PLACING PLIERS
- ELASTIC LIGATURE PLACING FORCEPS
- ELASTIC APPLICATORS & LIGATURE TUCKERS
- BRACKET POSITIONING INSTRUMENTS
- BRACKET & MEASURING INSTRUMENTS
- CROWN REMOVERS & FORCEPS FOR ARTICULATING PAPER
- CROWN SPREADERS
- CROWN INSTRUMENTS
- LABORATORY INSTRUMENTS



TM-116-1730
Brewer



TM-116-1731
Ivory



TM-116-1732
Ivory



TM-116-1733
Matrix Plier
17 cm



TM-116-1734
Stokes



TM-116-1735
Stokes / Palmer



TM-116-1736
Ivory

Rubber Dam Punches




Fig.1
TM-116-1737
Ainsworth



Fig.2
TM-116-1738
Ainsworth



Rubber Dam Instruments

RIGHT TEETH

TM-117-211 211 Winged TM-117-2A 2A Wingless TM-117-212 212 Wingless TM-117-29 29 Wingless

TM-117-51 51 Wingless TM-117-31 31 Wingless TM-117-30 30 Wingless TM-117-25 25 Wingless TM-117-W8A W8A Wingless TM-117-18 18 Wingless

TM-117-4 4 Winged TM-117-205 205 Winged TM-117-5 5 Winged TM-117-8A 8A Winged TM-117-201 201 Winged TM-117-14 14 Winged TM-117-51 51 Wingless TM-117-26 26 Wingless TM-117-W8A W8A Wingless TM-117-18 18 Wingless

TM-117-1 1 Winged TM-117-1A 1A Winged TM-117-207 207 Winged TM-117-208 208 Winged TM-117-2A 2A Wingless TM-117-22 22 Wingless TM-117-27 27 Wingless TM-117-W2A W2A Wingless

TM-117-11 11 Winged TM-117-10 10 Winged TM-117-3 3 Winged TM-117-203 203 Winged TM-117-204 204 Winged TM-117-202 202 Winged TM-117-201 201 Winged

TM-117-211 211 Winged TM-117-2 2 Wingless TM-117-22 22 Wingless

TM-117-211 211 Winged TM-117-210 210 Wingless TM-117-90N 90N Wingless

TM-117-6 6 Winged TM-117-9 9 Wingless

TM-117-211 211 Winged TM-117-212 212 Wingless

UPPER

PERMANENT

LOWER

LEFT TEETH

TM-117-208 208 Winged TM-117-207 207 Winged TM-117-206 206 Winged TM-117-22 22 Wingless TM-117-27 27 Wingless TM-117-29 29 Wingless


TM-117-W3 W3 Wingless TM-117-31 31 Wingless TM-117-30 30 Wingless TM-117-28 28 Wingless TM-117-25 25 Wingless TM-117-15 15 Wingless

TM-117-205 205 Winged TM-117-5 5 Winged TM-117-8A 8A Winged TM-117-201 201 Winged TM-117-14A 14A Winged TM-117-4 4 Wingless TM-117-8 8 Wingless TM-117-28 28 Wingless TM-117-24 24 Wingless

TM-117-207 207 Winged TM-117-209 209 Winged TM-117-0 0 Winged TM-117-00 00 Winged TM-117-2A 2A Wingless TM-117-1 1 Wingless TM-117-1A 1A Wingless


TM-117-200 200 Winged TM-117-202 202 Winged TM-117-14 14 Winged TM-117-56 56 Winged TM-117-7 7 Winged TM-117-7A 7A Winged TM-117-14A 14A Winged TM-117-56A 56A Winged

TM-117-29 29 Wingless TM-117-22 22 Wingless TM-117-27 27 Wingless TM-117-W2A W2A Wingless




TM-118-1741
Rubber Dam Clamps Tray


Rubber Dam Frames



TM-118-1742 13 cm
TM-118-1743 15 cm
Young



TM-118-1744
15 cm



TM-118-1745

Matrix Retainers



TM-118-1746
with cheek protector
5 mm



Fig. 8
TM-118-1747
Ivory
5 mm



Fig. 8
TM-118-1748
Ivory
6 mm



Fig. 8
TM-118-1749
Ivory
7 mm



TM-118-1750
Tofflemire junior



TM-118-1751
Tofflemire senior



TM-118-1752
Tofflemire universal



TM-118-1753
Nystrom
Left



TM-118-1754
Nystrom
Right



Fig. 12A
TM-118-1755
Ivory








Fig. 12
TM-118-1756
Ivory



TM-118-1757
Strip clamp
straight

Ligature & Separators Placing Pliers





TM-120-1764
Str 12 cm



TM-120-1765
Cvd 12 cm



TM-120-1766
14 cm



TM-120-1767
12 cm



TM-120-1768
14 cm



TM-120-1769
Elastomeric ligature
removing tweezer



Ligature plier
14 cm



Fig.1
TM-120-1770



Fig.2
TM-120-1771



Fig.3
TM-120-1772



Smaha
Ligature plier
14 cm



TM-120-1773
Str



TM-120-1774
Cvd



TM-120-1775
Ligature plier
14 cm



With tungsten
carbide insert



Elastic Applicators & Ligature Tuckers

Elastic Ligature Applicators

For easy placement of elastic ligatures under bracket tie wings



TM-121-1776



TM-121-1777

Elastomeric Ligature Removing Instruments

Multi-purpose instruments, ligature tucker & elastomeric ligature remover



TM-121-1778



TM-121-1779

Ligature Tuckers



TM-121-1780



TM-121-1781



TM-121-1782

Ligature Tuckers

Solid handle



TM-121-1783



TM-121-1784



TM-121-1785



Distal End Bending Instruments



TM-122-1786
For wire upto .22" (0.56 mm)

TM-122-1787
For wire upto .30" (0.76 mm)



TM-122-1788
For wire upto .22" (0.56 mm) Green

TM-122-1789
For wire upto .30" (0.76 mm) Green



Orange **TM-122-1790**

NiTi Ice Instrument 11.5 cm



TM-122-1791
Turns NiTi wire soft and flexible immediately for faster wire placement



Bracket Placer for Posterior Teeth



TM-122-1794

Bracket Positioning Gauges Anodized Aluminium



TM-122-1795 For .018" slot brackets



TM-122-1796 For .022" slot brackets



TM-122-1797
Boone
Universal

Adhesive Removers



TM-122-1792 **TM-122-1793**
Yellow

Bracket Placing Tweezers



TM-122-1798 13.5 cm



TM-122-1799 13.5 cm



TM-122-1800 11 cm

Bracket Tweezers with Slot Aligner



TM-122-1801 13.5 cm

Fig. 1



TM-122-1802 13.5 cm

Fig. 2



TM-122-1803 13.5 cm

Fig. 3



Bracket Measuring Instruments



TM-123-1804
Band Seater
15 cm



TM-123-1805
Mershon band pusher
14 cm



TM-123-1806
Band seating

Scaler & band pushing Instruments



TM-123-1807



TM-123-1808
Brown

Reverse action buccal tube tweezer




TM-123-1809 12.5 cm
TM-123-1810 13.5 cm

Buccal tube tweezer



TM-123-1811 12 cm
TM-123-1812 13.5 cm
Keat

Measuring Instruments




Bernklau
TM-123-1813 Orthodontic ruler
Perfect for all diagnostic measurements,
both intraoral and on the model



Ruler stainless steel
TM-123-1814 15 cm / 6"
TM-123-1815 30 cm / 12"




Beerendonk
TM-123-1816
Dental vernier with long & fine tips for intraoral measurements



TM-123-1817
Double sided, for inner & outer measurements



TM-123-1818
Movable Heads Bracket Positioning Gauge



Orthodontic bow compass
TM-123-1819

Turrets


For forming archwires, without torque,
for wires .016" to .022" (0.41mm to 0.56mm)




TM-123-1820
Made of Stainless steel



TM-123-1821
Made of anodized Aluminum
Disinfectable




TM-124-1822
Contents: 1 Crown Remover, Mounting Kit (Wrench & Key), 3 Tips




TM-124-1827
Contents: 1 Crown Remover, Mounting Kit (Wrench & Key), 3 Tips

Morrel Manual Crown Remover




TM-124-1832
Complete with 3 Tips

Morrel Manual Crown Remover with Spring




TM-124-1836
Complete with 3 Tips

Semi-Automatic Crown & Bridge Remover



TM-124-1840
Contents: 1 Crown Remover, Tightening Key, 3 Tips

Forceps For Articulating Paper



Crown Removers & Crown Spreaders

manual crown remover double ended



Fig. 1
TM-125-1849




Fig. 2
TM-125-1850



Fig. 3
TM-125-1851



TM-125-1848
Complete with 3 tips



TM-125-1852
Wynman
Crown Plier
with tungsten carbide inserts



TM-125-1853
Furrer Crown Plier
16 cm



TM-125-1854
Crown spreader



Fig. 3
TM-125-1855
Crown Spreader



Crown Spreader






Fig. 0
TM-125-1856



Fig. 1
TM-125-1857



Fig. 2
TM-125-1858



Furnace Tongs



TM-127-1867
31 cm



TM-126-1868
40 cm



TM-126-1869
56 cm



TM-126-1871 20 cm Str
TM-126-1872 20 cm Cvd
Plaster Cutter

Calipers



Fig.1
TM-126-1872
Boley
Measuring range
0-100 mm




Fig.2
TM-126-1873
Boley
Measuring range
0-100 mm



TM-126-1874
Iwanson for Wax
Measuring range
0-10 mm



TM-126-1875
Iwanson for Metal
Measuring range
0-10 mm



TM-126-1876
Micro Boley Gauge



Modeling saw & blades



TM-127-1877

Fret saw & blades



TM-127-1883

Saw blades for metal Pack of 12 pcs



TM-127-1878

0.25 x 3 mm Pack of 12 pcs



TM-127-1879

0.25 x 5 mm Pack of 12 pcs



TM-127-1880

0.12 x 5 mm Pack of 12 pcs



TM-127-1881

Soldering Tweezer
Str 16 cm



TM-127-1882

Soldering Tweezer
Cvd 16 cm



Fig.1

TM-127-1885
Riveting Mallet
22 cm



Fig.2

TM-127-1886
Riveting Mallet
22 cm



TM-127-1887

Anvil
Carbon Steel

Bib chains



Fig. 1
TM-129-1888
Ball chain



TM-129-1889
Plastic chain

X-Ray holders



Fig. 1
TM-129-1890



Fig. 6
TM-129-1891
6 Clips



Fig. 10
TM-129-1892
10 Clips




Fig. 14
TM-129-1893
14 Clips



TARIQ MAHMOOD
● ● ● ● ● enterprises

Page 130–131

o ARTICULATORS



TM-130-1894
Articulator ASA W/out Bearing Rod
Anatomical Articulator With Incisive
Pin, 4 Mounting Plates And Supporting
Back-pin. Available In Assorted Colors



TM-130-1895
ASA Articulator With Lateral Moment.
Made Of Aluminum Die Casting With
2 Mounting Plates. Also Available In
Vibratory Polished Process & Assorted colors



TM-130-1896
Keystone Deluxe Magnetic Denture
Articulator With Plates, Interchangeable,
Incisal Pin Included. Magnetic Split Casting
Ensures Easy & Exact Setting, Lateral
Excursions



TM-130-1897
Keystone Deluxe Universal Type-A Articulator,
Interchangeable, Included Incisal Pin, Lateral
Excursions



TM-130-1898
Anatomical Articulator Universal, Realized
For Easy Visibility To The Lingual Surface
& For A Perfect Stability Also In Reserved
Position. Also available in Assorted Colors



TM-130-1899
Chromium Plated Fork
Fork Crowned Bridge Nickel
Chromium Plated Articulator



TM-131-1900
Full Arch Hinged Articulator
Chromium Plated Hinge Articulator
For Mounting Full Arch Cases. Low
Back. Gypsum-free Articulation Enables
Clean & Accurate Mounting & Dismounting
In 30 Seconds



TM-131-1901
Ray Foster 57 Chrome Articulator
Chrome Brass Hinge Crowned Bridge
Articulator For Mounting Quadrant Cases.
Crown & Bridge Articulator Is Popular Ray
Foster # 57 Design, Eliminates The Need
To Break Off Models



TM-131-1902
Occlusor, Wire Articulator, Available
In 2 Sizes, Small & Large



TM-131-1903
Professional Occlusor, Vibratory
Polished Articulator.



TM-131-1904
Dental Articulator 2x-cresant Brand
Full Denture Chromium Plated Articulator.
Provides Additional Vertical Dimension For
Acceptance Of Higher Models. 49 Wide At
The Joint & Has A Vertical Span Between
Parallel Bows Of 25/8"



TM-131-1905
Hinge Articulator Full Arch
Plaine Line Nickle Chromium
Plated Articulator



TM-131-1906
Magic Art-I Articulator Lightweight



TM-131-1907
Mortimer-I Oculader Articulator



TM-131-1908
Mortimer-II Oculader Articulator

Page 132–149

- ORTHODONTIC PLIERS
- WIRE CUTTERS WITH TUNGSTEN CARBIDE INSERTS
- BONE CUTTERS
- INSTRUMENTS FOR "FIXED APPLIANCE TECHNIQUE"
- BONE RONGEURS



TM-132-1910
Flat nose wire bending plier
13 cm



TM-132-1911
Peeso
Round nose wire bending plier
13 cm



TM-132-1912
Peeso
Crown & collar contouring plier
13 cm



TM-132-1913
Aderer-Mini
11.5 cm



TM-132-1914
Aderer-Maxi
12 cm



TM-132-1915
Schwarz graduated
Optical style
13 cm



TM-132-1916
Schwarz
Optical style, round &
concave beaks 14 cm








TM-132-1917
Schwarz
Optical style, round &
concave beaks 14 cm




TM-132-1918
Mckellcps
Optical style, round &
concave beaks 16 cm





Orthodontic pliers





Replacement pads
TM-134-1932
TM-134-1933





Orthodontic pliers



TM-135-1940
Replacement tip stainless
440 for Omega or Tweed
loop forming pliers.
Tip diameters: 1.14mm, 1.52mm
& 1.91 mm



TM-135-1941
Omega 13 cm
Loop forming pliers with replaceable
stainless steel tip, three step tip has
1.14, 1.52 & 1.91 mm diameters




TM-135-1942
Tweed
Loop forming plier
13 cm



TM-135-1943
Angled



TM-135-1944
Straight



How
standard
4 mm wide tips
13 cm



TM-135-1945
Band Seating
Plier 13 cm



TM-135-1946
Nance
Loop forming plier
with four steps beak
13 cm



TM-135-1947
Tweed
Loop forming plier
13 cm








TM-135-1948
Weingart Mini
13 cm



TM-135-1949
Weingart
13 cm



TM-135-1950
Weingart Narrow
Narrow long tips for better
access to reach hard areas
13 cm





Orthodontic pliers



TM-137-1960
Adams Mini
12 cm



TM-137-1961
Adams Maxi
12.5 cm



TM-137-1962
Round nose plier
Fine pattern
13 cm



TM-137-1963
Adams Mini
13 cm



TM-137-1964
Adams Maxi
13 cm



TM-137-1965
Light wire bending
& loop forming Plier
13 cm



TM-137-1966
Angel
Bird beak plier
13 cm



TM-137-1967
Angel
Round nose wire
bending plier
13 cm



TM-137-1968
Angel
Round nose wire bending
& loop forming plier
13 cm






TM-138-1969
Serrated Laboratory
plier 14 cm



TM-138-1970
Smooth Laboratory
plier 14 cm



TM-138-1971
Laboratory plier
14.5 cm




TM-138-1972
Ribbon Arch plier
13 cm




TM-138-1973
Cone & Socket plier
13 cm



TM-138-1974
Goslee contouring
plier 13 cm



TM-138-1975
NiTi Cinch back
plier



TM-138-1976
Triplex Contouring
plier
14.5 cm



TM-138-1977
Wire bending
plier
14 cm



Orthodontic pliers



TM-139-1978
De La Rosa
Graduated
13 cm




TM-139-1979
De La Rosa
Smooth
13 cm



TM-139-1980
Wire bending plier
13 cm



TM-139-1981
Andresen wire
bending plier
13 cm



TM-139-1982
Hawley retainer
Plier
13 cm



TM-140-1983
Universal-Dolphin style
Wire bending & cutting plier
16 cm
● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")



TM-140-1984
Nance
Adams clasp forming plier
13 cm
Max. 0.70 mm (.028")



TM-140-1985
Universal-Taper style
Wire bending & cutting plier
16 cm
● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")



TM-140-1986
Schwarz
Arrow clasp forming plier
13 cm
Max. 0.70 mm (.028")



TM-140-1987
Waldasachs
Wire bending & cutting plier
15 cm
● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")



TM-140-1988
Schwarz
Arrow clasp forming plier
13 cm
Max. 0.70 mm (.028")



TM-140-1989
Universal
parallel jaw Plier
with side cutter
Max. 1.80 mm (.071")



TM-140-1990
Magnum
Universal facebow loop
forming & cutting plier
Max. 1.80 mm (.071")



Wire Cutters with Tungsten Carbide Inserts



TM-141-1991
Distal end cutter
with long handle
14 cm

●	Max. 0.30 mm to 0.50 mm (.012" to .020")
■	Max. 0.46 mm x 0.64 mm (.018" x .025")
⊗	Max. 0.36 mm to 0.53 mm (.0155" to .0215")
⊗	Max. 0.53 mm x 0.64 mm (.021" x .025")



TM-141-1992
Distal end cutter
with long handle
14 cm

●	Max. 0.30 mm to 0.50 mm (.012" to .020")
■	Max. 0.46 mm x 0.64 mm (.018" x .025")
⊗	Max. 0.36 mm to 0.53 mm (.0155" to .0215")
⊗	Max. 0.53 mm x 0.64 mm (.021" x .025")



TM-141-1993
Distal end cutter
14 cm

●	Max. 0.30 mm to 0.50 mm (.012" to .020")
■	Max. 0.46 mm x 0.64 mm (.018" x .025")
⊗	Max. 0.36 mm to 0.53 mm (.0155" to .0215")
⊗	Max. 0.53 mm x 0.64 mm (.021" x .025")



TM-141-1994
soft wire & ligature cutter
13 cm

.25 mm to .38 mm (.010" to .015")



TM-141-1995
soft wire & ligature cutter
13 cm

.25 mm to .38 mm (.010" to .015")



TM-141-1996
soft wire & ligature cutter
45° angled
13 cm

.25 mm to .38 mm (.010" to .015")

Heavy Box Joint Construction Cutters



TM-141-1997
Hard wire cutter, straight
14 cm

● Max. 0.30 mm to 0.50 mm (.012" to .020")

■ Max. 0.53 mm x 0.71 mm (.018" x .025")



TM-141-1998
Hard wire cutter, 15° angled
14 cm

● Max. 0.30 mm to 0.50 mm (.012" to .020")

■ Max. 0.53 mm x 0.71 mm (.022" x .028")



With Tungsten Carbide Inserts



TM-142-1999
How Standard
4 mm wide tips
13 cm



TM-142-2000
Angled

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-142-2001
Straight

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-142-2002
Weingart Mini
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-142-2003
Weingart
13 cm

● Max. 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-142-2004
Soft wire & ligature cutter
13 cm

□ .25 mm to .38 mm (.010" to .015")



TM-142-2005
Goldman-Fox
10 cm



Liston Bone cutter

TM-142-2006 14 cm Str M-142-2007 14 cm Cvd
TM-142-2008 17 cm Str M-142-2009 17 cm Cvd
TM-142-2010 19 cm Str M-142-2011 19 cm Cvd
TM-142-2012 22 cm Str M-142-2013 22 cm Cvd



Instruments for "Fixed Appliance Technique"

Standard distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2014
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Soft wire & ligature cutter, Mini

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2018
13 cm

25 mm to .38 mm (.010" to .015")

Mini head distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2015
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.64 mm (.018" x .025")
⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")

Soft wire and ligature cutter, standard

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2019
13 cm

25 mm to .38 mm (.010" to .015")

Distal end cutter with safety hold, mini head, long handles

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2016
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.64 mm (.018" x .025")
⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")

Soft wire & ligature cutter, 15° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2020
13 cm

25 mm to .38 mm (.010" to .015")

Flush-cut distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2017
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.64 mm (.018" x .025")

Soft wire & ligature cutter, 45° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-143-2021
13 cm

25 mm to .38 mm (.010" to .015")



Utility pliers



TM-144-2022
Weingart-Mini
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-144-2023
Weingart
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Utility pliers



TM-144-2024
How
Straight 13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



TM-144-2025
How
Angled 13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Wire bending & loop forming pliers



TM-144-2026
Angel (Bird beak)
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Light wire bending & loop forming pliers



TM-144-2027
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.41 mm x 0.56 mm (.016" x .022")

Light wire bending pliers with cutter Fitted with hard metal tips



TM-144-2028
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Light wire bending & loop forming pliers with cutter



TM-144-2029
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Light wire bending & loop forming pliers, grooved tip



TM-144-2030
Tweed
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")



Instruments for "Fixed Appliance Technique"

Hard wire cutter, straight

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-145-2031
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers angled



TM-145-2035
13 cm

Hard wire cutter, 15° angled

Non-autoclavable



TM-145-2032
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers straight



TM-145-2036
13 cm

Hard wire cutter, 15° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



TM-145-2033
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers with plastic tip



TM-145-2037
13 cm

Band removing pliers



TM-145-2034
13 cm

Adhesive removing pliers



TM-145-2038
13 cm

TM-145-2039 Replacement metal tip



Loop forming pliers with cutter



TM-146-2040

Tweed

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.53 mm x 0.64 mm (.021" x .025")

Archwire forming pliers, concave & convex tips



TM-146-2044

Hollow-Chop

13 cm

● 0.30 mm to 0.71 mm (.012" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Light wire bending pliers



TM-146-2041

Jarabak

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Multi-purpose pliers (bird beak, De la rosa & cutter in one)



TM-146-2045

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

3-jaw wire bending pliers



TM-146-2042

Aderer

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Distal end bending pliers



TM-146-2046

(Cinch back / tie back)

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Ligature forming pliers



TM-146-2043

13 cm

■ .25 mm to .38 mm (.010" to .015")

Distal end bending pliers



TM-146-2047

(Hammer head)

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



Instruments for "Fixed Appliance Technique"

Loop forming pliers



TM-147-2048

Nance

13 cm

● 0.36 mm to 0.71 mm (.014" to .028")

■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Torquing pliers set



0.18" 0.22"

TM-147-2052

13 cm

Detailing step pliers / Bayonet bending pliers - 1mm



TM-147-2049

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")

■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Rectangular arch forming & torquing pliers



TM-147-2053

Tweed / angel

13 cm

● 0.30 mm to 0.71 mm (.012" to .028")

■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Stop pliers V- Bend - 1mm



TM-147-2050

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")

■ Max. 0.43 mm x 0.56 mm (.017" x .022")

Surgical ball hook crimping pliers, angled



TM-147-2054

13 cm

Lingual / palatal arch forming pliers



TM-147-2051

13 cm





● 0.76 mm (.030") / 0.91 mm (.036")


Surgical ball hook crimping pliers, straight




TM-147-2055

13 cm





Friedmann-Mini
12.5 cm



TM-149-2071
Friedmann
14 cm



TM-149-2068

TM-149-2069

TM-149-2070



Friedmann-Micro
15 cm




Blumenthal
14.5 cm



TM-149-2072



TM-149-2073



TM-149-2074



TM-149-2075
Bayer
18 cm



TARIQ MAHMOOD
● ● ● ● ● enterprises


Page 150–154

o IMPRESSION TRAYS

Perforated without Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	TM-150-2100
L2	XS	60	49	TM-150-2101
L3	S	72	54	TM-150-2102
L4	M	74	57	TM-150-2103
L5	L	76	59	TM-150-2104
L6	XL	78	61	TM-150-2105



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	TM-150-2106
U2	XS	64	51	TM-150-2107
U3	S	66	53	TM-150-2108
U4	M	69	57	TM-150-2109
U5	L	71	60	TM-150-2110
U6	XL	73	62	TM-150-2111

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	TM-150-2112
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	TM-150-2113
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	TM-150-2114

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	TM-150-2115
L2	XS	60	49	TM-150-2116
L3	S	72	54	TM-150-2117
L4	M	74	57	TM-150-2118
L5	L	76	59	TM-150-2119
L6	XL	78	61	TM-150-2120



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	TM-150-2121
U2	XS	64	51	TM-150-2122
U3	S	66	53	TM-150-2123
U4	M	69	57	TM-150-2124
U5	L	71	60	TM-150-2125
U6	XL	73	62	TM-150-2126

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	TM-150-2127
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	TM-150-2128
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	TM-150-2129

Solid without Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	TM-150-2130
L2	XS	60	49	TM-150-2131
L3	S	72	54	TM-150-2132
L4	M	74	57	TM-150-2133
L5	L	76	59	TM-150-2134
L6	XL	78	61	TM-150-2135



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	TM-150-2136
U2	XS	64	51	TM-150-2137
U3	S	66	53	TM-150-2138
U4	M	69	57	TM-150-2139
U5	L	71	60	TM-150-2140
U6	XL	73	62	TM-150-2141

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	TM-150-2142
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	TM-150-2143
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	TM-150-2144

Solid with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	TM-150-2145
L2	XS	60	49	TM-150-2146
L3	S	72	54	TM-150-2147
L4	M	74	57	TM-150-2148
L5	L	76	59	TM-150-2149
L6	XL	78	61	TM-150-2150



Upper



No.	Size	A mm	B mm	Code
U1	XXS	61	47	TM-150-2151
U2	XS	64	51	TM-150-2152
U3	S	66	53	TM-150-2153
U4	M	69	57	TM-150-2154
U5	L	71	60	TM-150-2155
U6	XL	73	62	TM-150-2156

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	TM-150-2157
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	TM-150-2158
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	TM-150-2159




Impression Trays Regular

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XL	80	63	TM-151-2160
L2	L	78	59	TM-151-2161
L3	M	76	57	TM-151-2162
L4	NM	72	57	TM-151-2163
L5	S	69	50	TM-151-2164



Upper




No.	Size	A mm	B mm	Code
U1	XL	80	63	TM-151-2165
U2	L	78	59	TM-151-2166
U3	M	76	57	TM-151-2167
U4	NM	72	57	TM-151-2168
U5	S	69	50	TM-151-2169

Set	Contents	Code
10 Pcs	L1,L2,L3,L4,L5, U1,U2,U3,U4,U5	TM-151-2170

Solid with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XL	80	63	TM-151-2171
L2	L	78	59	TM-151-2172
L3	M	76	57	TM-151-2173
L4	NM	72	57	TM-151-2174
L5	S	69	50	TM-151-2175



Upper



No.	Size	A mm	B mm	Code
U1	XL	80	63	TM-151-2176
U2	L	78	59	TM-151-2177
U3	M	76	57	TM-151-2178
U4	NM	72	57	TM-151-2179
U5	S	69	50	TM-151-2180


Set	Contents	Code
10 Pcs	L1,L2,L3,L4,L5, U1,U2,U3,U4,U5	TM-151-2181

Impression Trays For Children

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	S	54	47	TM-151-2182
L2	M	62	51	TM-151-2183
L3	L	70	55	TM-151-2184



Upper




No.	Size	A mm	B mm	Code
U1	S	55	43	TM-151-2185
U2	M	60	45	TM-151-2186
U3	L	63	49	TM-151-2187

Set	Contents	Code
6 Pcs	L1,L2,L3,U1,U2,U3	TM-151-2188



Lower


Solid with Rim



No.	Size	A mm	B mm	Code
L1	S	54	47	TM-151-2189
L2	M	62	51	TM-151-2190
L3	L	70	55	TM-151-2191



Upper




No.	Size	A mm	B mm	Code
U1	S	55	43	TM-151-2192
U2	M	60	45	TM-151-2193
U3	L	63	49	TM-151-2194

Set	Contents	Code
6 Pcs	L1,L2,L3,U1,U2,U3	TM-151-2195

Impression Trays Regular Ehrlicke Pattern Perforated




Lower



No.	Size	A mm	B mm	Code
UB0	XS	65	46	TM-151-2196
UB1	S	67	50	TM-151-2197
UB2	M	76	54	TM-151-2198
UB3	L	79	56	TM-151-2199
UB4	XL	79	61	TM-151-2200



Upper



No.	Size	A mm	B mm	Code
UB0	XS	59	52	TM-151-2201
UB1	S	66	52	TM-151-2202
UB2	M	71	62	TM-151-2203
UB3	L	71	64	TM-151-2204
UB4	XL	80	70	TM-151-2205

Set	Contents	Code
10 Pcs	LB0,LB1,LB2,LB3,LB4, UB0,UB1,UB2,UB3,UB4	TM-151-2206



Impression Trays

Impression Trays Edentulous

Perforated with Rim



Lower

No.	Size	A mm	B mm	Code
L1	XL	84	65	TM-152-2207
L2	L	77	64	TM-152-2208
L3	M	76	60	TM-152-2209
L4	S	68	57	TM-152-2210



Upper

				
No.	Size	A mm	B mm	Code
U1	XL	80	62	TM-152-2211
U2	L	77	60	TM-152-2212
U3	M	70	60	TM-152-2213
U4	S	72	56	TM-152-2214
Set	Contents			Code
8 Pcs	L1,L2,L3,L4,U1,U2,U3,U4			TM-152-2215

Solid with Rim



Lower

No.	Size	A mm	B mm	Code
L1	XL	84	65	TM-152-2216
L2	L	77	64	TM-152-2217
L3	M	76	60	TM-152-2218
L4	S	68	57	TM-152-2219



Upper

				
No.	Size	A mm	B mm	Code
U1	XL	80	62	TM-152-2220
U2	L	77	60	TM-152-2221
U3	M	70	60	TM-152-2222
U4	S	72	56	TM-152-2223
Set	Contents			Code
8 Pcs	L1, L2, L3, L4, U1, U2, U3, U4			TM-152-2224

Impression Trays Anterior Depressed

Perforated with Rim



Lower

No.	Size	A mm	B mm	Code
L1	S	65	48	TM-152-2225
L2	M	72	57	TM-152-2226
L3	L	75	57	TM-152-2227



Upper

				
No.	Size	A mm	B mm	Code
U1	S	64	56	TM-152-2228
U2	M	67	60	TM-152-2229
U3	L	72	63	TM-152-2230
Set	Contents			Code
6 Pcs	L1,L2,L3,U1,U2,U3			TM-152-2231

Solid with Rim



Lower

No.	Size	A mm	B mm	Code
L1	S	65	48	TM-152-2232
L2	M	72	57	TM-152-2233
L3	L	75	57	TM-152-2234



Upper

				
No.	Size	A mm	B mm	Code
U1	S	64	56	TM-152-2235
U2	M	67	60	TM-152-2236
U3	L	72	63	TM-152-2237
Set	Contents			Code
6 Pcs	L1,L2,L3,U1,U2,U3			TM-152-2238



Perforated with Rim



Fig. 30
TM-153-2239
Upper Left / Lower Right



Fig. 31
TM-153-2240
Lower Left / Upper Right



Fig. 32
TM-153-2241
Lower Jaw



Fig. 33
TM-153-2242
Upper Jaw



Fig. 99
TM-153-2243
Lateral

Solid with Rim



Fig. 30
TM-153-2244
Upper Left / Lower Right



Fig. 31
TM-153-2245
Lower Left / Upper Right



Fig. 32
TM-153-2246
Lower Jaw



Fig. 33
TM-153-2247
Upper Jaw



Fig. 99
TM-153-2248
Lateral

Perforated without Rim



Lower Left / Upper Right		
No.	Size	Code
Fig. L2D	S	TM-153-2249
Fig. L4D	M	TM-153-2250
Fig. L6D	L	TM-153-2251



Upper Left / Lower Right		
No.	Size	Code
Fig. L2S	S	TM-153-2255
Fig. L4S	M	TM-153-2256
Fig. L6S	L	TM-153-2257

Perforated with Rim



Lower Left / Upper Right		
No.	Size	Code
Fig. L2D	S	TM-153-2252
Fig. L4D	M	TM-153-2253
Fig. L6D	L	TM-153-2254



Upper Left / Lower Right		
No.	Size	Code
Fig. L2S	S	TM-153-2258
Fig. L4S	M	TM-153-2259
Fig. L6S	L	TM-153-2260

Carvers for silicone impression materials




TM-154-2261




TM-154-2262


Gingival cord packers




TM-154-2263 Serrated



TM-154-2265 Serrated



TM-154-2264 Smooth



TM-154-2266 Serrated



TARIQ MAHMOOD
enterprises

Page 155–164

- SPIRIT LAMPS & JARS
- BOWLS & KIDNEY TRAYS
- INSTRUMENTS TRAYS
- INSTRUMENTS BOXES
- STERILIZING BOXES & INSTRUMENTS TRAY
- MEASURING JUGS, BUCKETS & IRRIGATORS
- NEEDLE CASES, PETRI DISHES & FUNNELS
- INSTRUMENTS TRAYS & ACCESSORIES FOR BURS
- ENDODONTIC TRAYS & INSTRUMENTS



TM-155-2300 Spirit Lamp 120 cc



TM-155-2301 Spirit Lamp 120 cc



TM-155-2302 Spirit Lamp 120 cc



TM-155-2303 Spare wicks Flat 12mm (pack of 5)



TM-155-2304 Spare wicks Round Ø4mm (pack of 5)



Forceps Jars

	dia.	height
TM-155-2305	30	x 90 mm
TM-155-2306	50	x 100 mm
TM-155-2307	50	x 130 mm
TM-155-2308	50	x 175 mm
TM-155-2309	75	x 100 mm
TM-155-2310	75	x 130 mm
TM-155-2311	75	x 175 mm
TM-155-2312	100	x 175 mm
TM-155-2313	100	x 200 mm



Dressing Jars With Lid & Knob

	dia.	height
TM-155-2314	75	mm x 60 mm
TM-155-2315	75	mm x 75 mm
TM-155-2316	100	mm x 60 mm
TM-155-2317	100	mm x 100 mm
TM-155-2318	125	mm x 80 mm
TM-155-2319	125	mm x 125 mm
TM-155-2320	150	mm x 150 mm



Dressing Jars / Cotton Containers

	dia.	height
TM-155-2321	80	x 60 mm
TM-155-2322	80	x 80 mm
TM-155-2323	100	x 60 mm
TM-155-2324	100	x 80 mm
TM-155-2325	100	x 100 mm
TM-155-2326	120	x 80 mm
TM-155-2327	120	x 100 mm
TM-155-2328	120	x 120 mm
TM-155-2329	150	x 100 mm
TM-155-2330	150	x 150 mm
TM-155-2331	180	x 150 mm
TM-155-2332	180	x 180 mm
TM-155-2333	200	x 200 mm
TM-155-2334	240	x 240 mm



Gallipot	
	vol.
TM-156-2335	50 ml
TM-156-2336	100 ml
TM-156-2337	150 ml
TM-156-2338	225 ml



Stainless Steel Round Bowls			
	dia.	height	vol.
TM-156-2339	Ø 40 mm x 19 mm		0.05L
TM-156-2340	Ø 60 mm x 30 mm		0.07L
TM-156-2341	Ø 80 mm x 40 mm		0.14L
TM-156-2342	Ø 100 mm x 50 mm		0.35L
TM-156-2343	Ø 120 mm x 60 mm		0.45L
TM-156-2344	Ø 150 mm x 70 mm		0.70L
TM-156-2345	Ø 180 mm x 90 mm		1.10L
TM-156-2346	Ø 200 mm x 100 mm		1.60L



Stainless Steel Bowls			
	dia.	height	vol.
TM-156-2347	Ø 40 mm x 19 mm		0.05L
TM-156-2348	Ø 60 mm x 30 mm		0.07L
TM-156-2349	Ø 80 mm x 40 mm		0.14L
TM-156-2350	Ø 100 mm x 50 mm		0.35L
TM-156-2351	Ø 120 mm x 60 mm		0.45L
TM-156-2352	Ø 150 mm x 70 mm		0.70L
TM-156-2353	Ø 180 mm x 90 mm		1.10L
TM-156-2354	Ø 200 mm x 100 mm		1.60L



Wash Basins	
	dia.
TM-156-2355	300 mm
TM-156-2356	320 mm
TM-156-2357	350 mm
TM-156-2358	400 mm



Medicine Cup		
	dia.	height
TM-156-2359	80 x 35 mm	
TM-156-2360	125 x 60 mm	



Medicine Cup (Graduated)		
	Size (mm)	Cap.(ml)
TM-156-2361	2 Oz.	60 ml



Kidney Trays With LID	
	dia.
TM-156-2362	150 mm
TM-156-2363	200 mm
TM-156-2364	250 mm
TM-156-2365	300 mm



Kidney Trays	
	dia.
TM-156-2366	150 mm
TM-156-2367	200 mm
TM-156-2368	250 mm
TM-156-2369	300 mm



Kidney Trays Shallow	
	dia.
TM-156-2370	150 mm
TM-156-2371	200 mm
TM-156-2372	250 mm
TM-156-2373	300 mm



Scalers Trays	
TM-157-2374	190 x 100 x 17 mm
TM-157-2375	200 x 100 x 25 mm



Scalers Trays	
TM-157-2376	185 x 80 x 15 mm



Instruments Trays With LID	
TM-157-2377	150 x 100 x 50 mm
TM-157-2378	200 x 100 x 50 mm
TM-157-2379	220 x 120 x 50 mm
TM-157-2380	240 x 140 x 50 mm
TM-157-2381	200 x 150 x 50 mm
TM-157-2382	250 x 200 x 50 mm
TM-157-2383	300 x 250 x 50 mm
TM-157-2384	300 x 200 x 50 mm
TM-157-2385	350 x 250 x 50 mm
TM-157-2386	400 x 300 x 50 mm
TM-157-2387	450 x 300 x 50 mm
TM-157-2388	450 x 350 x 50 mm



Instruments Trays Without LID	
TM-157-2389	150 x 100 x 50 mm
TM-157-2390	200 x 100 x 50 mm
TM-157-2391	220 x 120 x 50 mm
TM-157-2392	240 x 140 x 50 mm
TM-157-2393	200 x 150 x 50 mm
TM-157-2394	250 x 200 x 50 mm
TM-157-2395	300 x 250 x 50 mm
TM-157-2396	300 x 200 x 50 mm
TM-157-2397	350 x 250 x 50 mm
TM-157-2398	400 x 300 x 50 mm
TM-157-2399	450 x 300 x 50 mm
TM-157-2400	450 x 350 x 50 mm



Instruments Trays	
TM-158-2401	210 x 160 x 10 mm
TM-158-2402	240 x 180 x 10 mm
TM-158-2403	310 x 150 x 10 mm
TM-158-2404	310 x 210 x 10 mm
TM-158-2405	350 x 240 x 10 mm
TM-158-2406	400 x 270 x 10 mm
TM-158-2407	420 x 280 x 10 mm



Mayo Trays Perforated	
TM-158-2408	150 x 100 x 50 mm
TM-158-2409	200 x 100 x 50 mm
TM-158-2410	220 x 120 x 50 mm
TM-158-2411	240 x 140 x 50 mm
TM-158-2412	200 x 150 x 50 mm
TM-158-2413	250 x 200 x 50 mm
TM-158-2414	300 x 200 x 50 mm
TM-158-2415	300 x 250 x 50 mm
TM-158-2416	350 x 250 x 50 mm
TM-158-2417	400 x 300 x 50 mm
TM-158-2418	450 x 300 x 50 mm
TM-158-2419	450 x 350 x 50 mm



DIN Baskets perforated sides	
TM-158-2420	240 x 240 x 50 mm
TM-158-2421	240 x 240 x 100 mm
TM-158-2422	480 x 260 x 50 mm
TM-158-2423	520 x 260 x 80 mm
TM-158-2424	520 x 260 x 100 mm



Wiremesh Baskets perforated	
TM-158-2425	240 x 240 x 50 mm
TM-158-2426	260 x 240 x 100 mm
TM-158-2427	480 x 260 x 50 mm
TM-158-2428	520 x 260 x 80 mm
TM-158-2429	520 x 260 x 100 mm



Instrument Boxes With Side Handle	
TM-159-2430	150 x 100 x 50 mm
TM-159-2431	200 x 100 x 50 mm
TM-159-2432	220 x 120 x 50 mm
TM-159-2433	240 x 140 x 50 mm
TM-159-2434	200 x 150 x 50 mm
TM-159-2435	250 x 200 x 50 mm
TM-159-2436	300 x 200 x 50 mm
TM-159-2437	300 x 250 x 50 mm
TM-159-2438	400 x 300 x 50 mm
TM-159-2439	450 x 300 x 50 mm
TM-159-2440	450 x 350 x 50 mm



Instrument Boxes With LID	
TM-159-2441	150 x 100 x 50 mm
TM-159-2442	200 x 100 x 50 mm
TM-159-2443	220 x 120 x 50 mm
TM-159-2444	240 x 140 x 50 mm
TM-159-2445	200 x 150 x 50 mm
TM-159-2446	250 x 200 x 50 mm
TM-159-2447	300 x 200 x 50 mm
TM-159-2448	300 x 250 x 50 mm
TM-159-2449	400 x 300 x 50 mm
TM-159-2450	400 x 300 x 50 mm
TM-159-2451	450 x 350 x 50 mm



Perforated Instruments Boxes	
TM-159-2452	170 x 070 x 35 mm
TM-159-2453	180 x 080 x 40 mm
TM-159-2454	200 x 100 x 40 mm
TM-159-2455	200 x 100 x 50 mm
TM-159-2456	200 x 100 x 60 mm
TM-159-2457	220 x 120 x 60 mm
TM-159-2458	250 x 120 x 60 mm
TM-159-2459	300 x 150 x 60 mm
TM-159-2460	300 x 200 x 80 mm
TM-159-2461	400 x 160 x 80 mm
TM-159-2462	400 x 200 x 100 mm
TM-159-2463	500 x 200 x 120 mm



Perforated Instruments Boxes	
TM-159-2464	170 x 070 x 35 mm
TM-159-2465	180 x 080 x 40 mm
TM-159-2466	200 x 100 x 40 mm
TM-159-2467	200 x 100 x 50 mm
TM-159-2468	200 x 100 x 60 mm
TM-159-2469	220 x 120 x 60 mm
TM-159-2470	250 x 120 x 60 mm
TM-159-2471	300 x 150 x 60 mm
TM-159-2472	300 x 200 x 80 mm
TM-159-2473	400 x 160 x 80 mm
TM-159-2474	400 x 200 x 100 mm
TM-159-2475	500 x 200 x 120 mm



Catheter trays with cover

TM-160-2476	100 x 200 x 50 mm
TM-160-2477	100 x 300 x 50 mm
TM-160-2478	120 x 220 x 50 mm



Instruments Boxes Perforated

TM-160-2479	180 x 80 x 30 mm
TM-160-2480	180 x 80 x 40 mm
TM-160-2481	180 x 80 x 50 mm
TM-160-2482	200 x 90 x 30 mm
TM-160-2483	200 x 90 x 40 mm
TM-160-2484	200 x 90 x 50 mm
TM-160-2485	250 x 100 x 40 mm
TM-160-2486	250 x 100 x 50 mm
TM-160-2487	250 x 100 x 60 mm



Thermometer Jar


	dia. height
TM-160-2488	Ø 25 x 100 mm



Instruments Boxes With Filters

TM-160-2492	170 x 070 x 35 mm
TM-160-2493	180 x 080 x 40 mm
TM-160-2494	200 x 100 x 40 mm
TM-160-2495	200 x 100 x 50 mm
TM-160-2496	200 x 100 x 60 mm
TM-160-2497	220 x 120 x 60 mm
TM-160-2498	250 x 120 x 60 mm
TM-160-2499	300 x 150 x 60 mm
TM-160-2500	300 x 200 x 80 mm
TM-160-2501	400 x 160 x 80 mm
TM-160-2502	400 x 200 x 100 mm
TM-160-2503	500 x 200 x 120 mm

Sterilizing Forceps



Cheatele Forceps

TM-160-2489 27 cm

Collin Forceps

TM-160-2490 20 cm
TM-160-2491 30 cm



Measuring Jugs , Buckets Irrigators



Irrigators	
	vol.
TM-161-2504	1.0L
TM-161-2505	1.5L
TM-161-2506	2.0L



Measuring Jugs	
	vol.
TM-161-2507	250 ml
TM-161-2508	500 ml
TM-161-2509	1000 ml
TM-161-2510	2000 ml



Measuring Jugs	
	vol.
TM-161-2511	250 ml
TM-161-2512	500 ml
TM-161-2513	1000 ml
TM-161-2514	2000 ml



Buckets With LID		Buckets Without LID	
	vol.		vol.
TM-161-2515	8.0L	TM-161-2518	8.0L
TM-161-2516	10.0L	TM-161-2519	10.0L
TM-161-2517	15.0L	TM-161-2520	0L



Sputum Cup	
	dia. height
TM-161-2521	Ø 75 x 80 mm
TM-161-2522	Ø 100 x 50 mm
TM-161-2523	Ø 100 x 80 mm



Ointment Jar		
	dia. height	Vol.
TM-161-2524	Ø 75 x 65 mm	250ml



Needle Case

TM-162-2525 210 x 125 x 17 mm



Needle Case Sauerbruch

TM-162-2526 210 x 125 x 17 mm



Needle Case

TM-162-2527 150 x 90 x 10 mm



Needle Cases Perforated

TM-162-2528 50 x 25 x 4.5 mm

TM-162-2529 65 x 42 x 8 mm

TM-162-2530 75 x 25 x 4.5 mm



Needle Cases Plain

TM-162-2531 50 x 25 x 4.5 mm

TM-162-2532 65 x 42 x 8 mm

TM-162-2533 75 x 25 x 4.5 mm



TM-162-2534 Needle Case Complete



Funnels

dia.

TM-162-2535 100 mm

TM-162-2536 120 mm

TM-162-2537 150 mm

TM-162-2538 200 mm

TM-162-2539 240 mm



Petri Dishes

dia. height

TM-162-2540 60 x 20 mm

TM-162-2541 90 x 15 mm

TM-162-2542 100 x 20 mm

TM-162-2543 120 x 20 mm



Instruments Trays & Accessories for Burs



Chico

Stainless Steel for 5 instruments

TM-163-2544 180 x 60 x 22 mm



Medio

Stainless Steel for 10 instruments

TM-163-2545 180 x 125 x 22 mm



Maxi

Stainless Steel for 20 instruments

TM-163-2546 280 x 180 x 25 mm



Maxi Set 1

Stainless Steel rack for 15 instruments

TM-163-2547 288 x 187 x 40 mm



Maxi Set 2

Stainless Steel for 15 instruments

TM-163-2548 288 x 187 x 40 mm



Maxi Set 3

Stainless Steel for 15 instruments

TM-163-2549 288 x 187 x 40 mm



Medio

Stainless Steel Perforated Tray with cover

TM-163-2550 187 x 145 x 25 mm



Medio

Stainless Steel Solid Tray with cover

TM-163-2551 187 x 145 x 25 mm



Burs Disinfecting Box Stainless Steel

TM-163-2552 90 mm x 40 mm



Mesh Burs Holder

TM-163-2553 35 mm x 150 mm



TM-163-2554 Bur Holding Tweezer 16.5 cm



Perforated Instruments Boxes

TM-164-2555	170 x 070 x 35 mm
TM-164-2556	180 x 080 x 40 mm
TM-164-2557	200 x 100 x 40 mm
TM-164-2558	200 x 100 x 50 mm
TM-164-2559	200 x 100 x 60 mm
TM-164-2560	220 x 120 x 60 mm
TM-164-2561	250 x 120 x 60 mm
TM-164-2562	300 x 150 x 60 mm
TM-164-2563	300 x 200 x 80 mm
TM-164-2564	400 x 160 x 80 mm
TM-164-2565	400 x 200 x 100 mm
TM-164-2566	500 x 200 x 120 mm



Stainless Steel endodontic tray with cover

TM-164-2567 250 x 115 x 15 mm



Instrument Box perforated

TM-164-2568 180 x 90 x 30 mm



TM-164-2569 Endometer Laser etched




TM-164-2570 Lester Pumice Cup



TM-164-2571 Endo Box Round



TM-164-2572 Endo Box



TM-164-2573 Endo Box Trio



TM-164-2574 Cotton Holder with Lid



TM-164-2575 Cotton Holder without Lid



TM-164-2576 Cotton Holder Estrella without Lid

Page 165–170

- ORAL SURGICAL BURS
- TRIMMER BURS
- LINDEMANN SIDE CUTTING BURS
- ENT BURS





Toller Round Head Bur

Tungsten Carbide
with Blood Ring



1/1



Length	50mm	70mm
Head size	2.3mm	2.3mm
Code No.	TM-165-2700	TM-165-2701



Sterile Packed in boxes of 10

Toller Taper Fissure Bur

Tungsten Carbide cross cut
with Blood Ring



1/1



Length	50mm	70mm
Head size	1.6mm	1.6mm
Code No.	TM-165-2702	TM-165-2703

Sterile Packed in boxes of 10

Apicectomy Bur

44mm length



1/1



Head size	0.8mm
Code No.	TM-165-2704

Sterile Packed in boxes of 10

Micro Head Apicectomy Bur

10mm length
Latch Fitting



1/1



Head size	1.0mm
Code No.	TM-165-2705







Sterile Packed in boxes of 5

Round Head Bur








Plain cut
44mm length



1/1



Head size	1.0mm	1.2mm	1.4mm	1.6mm	1.8mm	2.1mm
Code No.	TM-165-2706	TM-165-2707	TM-165-2708	TM-165-2709	TM-165-2710	TM-165-2711



2.3mm	2.5mm	2.7mm	3.1mm	4.0mm	5.0mm	6.0mm
TM-165-2712	TM-165-2713	TM-165-2714	TM-165-2715	TM-165-2716	TM-165-2717	TM-165-2718

Sterile Packed in boxes of 10



TARIQ MAHMOOD
enterprises


Oral Surgery Burs

Flat Fissure Bur

Cross cut
44mm length



1/1



Head size 1mm 1.4mm 1.6mm 1.8mm 2.1mm 2.3mm

Code No. **TM-166-2719** **TM-166-2720** **TM-166-2721** **TM-166-2722** **TM-166-2723** **TM-166-2724**


Sterile Packed in boxes of 10

Tapered Fissure Bur

Cross cut
44mm length



1/1



Head size 1mm 1.4mm 1.6mm 1.8mm 2.1mm 2.3mm

Code No. **TM-166-2725** **TM-166-2726** **TM-166-2727** **TM-166-2728** **TM-166-2729** **TM-166-2730**

Sterile Packed in boxes of 10

Tapered Fissure Bur

Cross cut
50mm length



Head size 1.6mm

Code No. **TM-166-2731**

Sterile Packed in boxes of 10




Bud Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

TM-167-2732 TM-167-2733 TM-167-2734

70mm length Code No.

TM-167-2735 TM-167-2736 TM-167-2737

44mm Sterile Packed in boxes of 10

70mm Sterile Packed in boxes of 5


Cylinder Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

TM-167-2738 TM-167-2739 TM-167-2740

70mm length Code No.

TM-167-2741 TM-167-2742 TM-167-2743

44mm Sterile Packed in boxes of 10

70mm Sterile Packed in boxes of 5


Pear Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

TM-167-2744 TM-167-2745 TM-167-2746

70mm length Code No.

TM-167-2747 TM-167-2748 TM-167-2749

44mm Sterile Packed in boxes of 10


70mm Sterile Packed in boxes of 5



TARIQ MAHMOOD
enterprises

Lindemann Side Cutting Burs


Lindemann Side Cutting Burs



Sterile Packed in boxes of 5

Bone Trephines

Bone Trephines



Sterile Packed in boxes of 5

Sterile Packed in boxes of 5




ENT Cutting Bur

70mm length




1/1



Head size 0.6mm 0.7mm 0.8mm 1mm 1.4mm 1.8mm 2.3mm

Code No. **TM-169-2757** **TM-169-2758** **TM-169-2759** **TM-169-2760** **TM-169-2761** **TM-169-2762** **TM-169-2763**



2.7mm 3.1mm 3.5mm 4mm 4.5mm 5mm 6mm 7mm

TM-169-2764 **TM-169-2765** **TM-169-2766** **TM-169-2767** **TM-169-2768** **TM-169-2769** **TM-169-2770** **TM-169-2771**

Sterile Packed in boxes of 5


ENT Allport Cutting Bur

70mm length

6 Blades




1/1



Head size 0.6mm 0.7mm 0.8mm 1mm 1.4mm 1.8mm 2.3mm

Code No. **TM-169-2772** **TM-169-2773** **TM-169-2774** **TM-169-2775** **TM-169-2776** **TM-169-2777** **TM-169-2778**



2.7mm 3.1mm 3.5mm 4mm 4.5mm 5mm 6mm 7mm

TM-169-2779 **TM-169-2780** **TM-169-2781** **TM-169-2782** **TM-169-2783** **TM-169-2784** **TM-169-2785** **TM-169-2786**

Sterile Packed in boxes of 5



TARIQ MAHMOOD
enterprises


ENT Burs

ENT Polishing Bur

70mm length




1/1



Head size 0.6mm 0.7mm 0.8mm 1mm 1.4mm 1.8mm 2.3mm

Code No. **TM-170-2787** **TM-170-2788** **TM-170-2789** **TM-170-2790** **TM-170-2791** **TM-170-2792** **TM-170-2793**



2.7mm 3.1mm 3.5mm 4mm 4.5mm 5mm 6mm 7mm

TM-170-2794 **TM-170-2795** **TM-170-2796** **TM-170-2797** **TM-170-2798** **TM-170-2799** **TM-170-2800** **TM-170-2801**


Sterile Packed in boxes of 5

ENT Diamond Bur

70mm length




1/1



Head size 0.6mm 0.7mm 0.8mm 1mm 1.4mm 1.8mm 2.3mm

Code No. **TM-170-2802** **TM-170-2803** **TM-170-2804** **TM-170-2805** **TM-170-2806** **TM-170-2807** **TM-170-2808**



2.7mm 3.1mm 3.5mm 4mm 4.5mm 5mm 6mm 7mm

TM-170-2809 **TM-170-2810** **TM-170-2811** **TM-170-2812** **TM-170-2813** **TM-170-2814** **TM-170-2815** **TM-170-2816**

Sterile Packed in boxes of 5

Registered Office:



77, Aziz Shaheed Colony
Sialkot Cantt
51300-Pakistan

M: +92 (0) 324 8472499

E: info@tmentp.com

Web: www.tmentp.com

Company Registration No. A-22629
National Tax Number (NTN) 4174224



ISO 9001:2015
ISO 13485:2016